



FUTURE

MOBILITY HEALTHCARE INC.

**POWER COMPONENTS
FOR POWER
WHEELCHAIR**

PARTS MANUAL

FUTURE MOBILITY HEALTHCARE APPRECIATES YOUR BUSINESS.
OUR CUSTOMER SERVICE REPRESENTATIVES CAN BE CONTACTED AT:

**TO PLACE ORDERS:
(905) 671-1661**

FOR YOUR CONVENIENCE YOU MAY ALSO PLACE ORDERS BY FAX AT:

FAX: (905) 671-3377

PLEASE VISIT OUR WEBSITE AT:

www.futuremobility.ca

*Part prices and configurations included in this Parts Manual are subject to change without prior notice.


* All prices are shown as per unit individual price.

Contents

POWER	5
POWER COMPONENTS AND ACCESSORIES	5

POWER

POWER COMPONENTS AND ACCESSORIES



POWER COMPONENTS AND ACCESSORIES - TABLE

ITEM	QTY	NOTES	PART NO.	DESCRIPTION	PRICE/ UNIT
1	1	A,E	17001296	PG VR-2 POWER MODULE PM 60 /D50681	\$576.00
		A,F	13001414	PG VR-2 POWER MODULE PM 90 /D51080	\$800.00
	1	B,E	17001411	PG R-NET POWER MODULE PM 80-1.2M /D50615	\$1,395.00
		B,F	1026475	PG R-NET POWER MODULE PM 120-1.2M /D50600	\$1,550.00
2	1	A	1011007	VR-2 JOYSTICK MODULE-JSM-CN /17001295	\$468.00
	1	B	17001412	R-NET JOYSTICK MODULE-JSM-SW /D50616	\$930.00
3	1	C	1015457	22 NF GEL BATTERY	\$350.00
		C	1026474	24 NF GEL BATTERY	\$450.00
4	1	D	17001365	VR-2 BATTERY CONNECTOR KIT	\$10.00
	1	D	17001366	VR-2 MOTOR CONNECTOR KIT	\$10.00
	1	D	17001415	R-NET BATTERY CONNECTOR KIT	\$10.00
	1	D	17001414	R-NET MOTOR CONNECTOR KIT	\$22.00
5	1	E	16001182	2-POLE MOTOR UNIT, RT	\$798.00
	1	E	16001181	2-POLE MOTOR UNIT, LT	\$798.00
	1	F	16001400	4-POLE MOTOR UNIT, RT	\$950.00
	1	F	16001399	4-POLE MOTOR UNIT, LT	\$950.00
6	1	E	17001361	SONEIL BATTERY CHARGER W POWER CORD	\$294.00
		F	17001423	BATTERY CHARGER-6AMP ONBOARD(2412-SRF)	\$394.00
7	1		17001406	VR-2 PG PP1A PROGRAMMER	\$875.00
	1		1026477	R-NET PG DONGLE PROGRAMMER	\$500.00

NOTES:

- A- VR-2 POWER MODULE AND VR-2 JOYSTICK ARE USED TOGETHER**
- B- R-NET POWER MODULE AND R-NET JOYSTICK ARE USED TOGETHER**
- C- PRICE FOR ONE UNIT ONLY**
- D- SPECIFY THE TYPE OF POWER CHAIR FOR THE CORRESPONDING CONNECTOR**
- E- FOR MAGELLAN POWER CHAIR**
- F- FOR VIKING POWER CHAIR**