



<b>Account Information</b>		<b>Order Date:</b> _____	
Account #:	_____	P.O.#	_____
Business Name:	_____	Phone:	_____
Ship to Address:	_____	City:	_____
Fax:	_____	Province:	_____
Contact:	_____	Postal Code:	_____
<b>TAG FOR:</b> _____			

**① FRAME TYPE**

- R  Orion II - WMA5FM005----- 101 **\$2,700**
- M  Power Tilt ----- 101MA31 **\$1,600**

**② FRAME FINISH**

- Black ----- 101C01 **No Charge**
- Silver ----- 101C02 **No Charge**
- Custom Colours**
- Blue----- 101C03 **\$150**
- Green----- 101C04 **\$150**
- Orange----- 101C05 **\$150**
- Purple----- 101C06 **\$150**
- Red----- 101C07 **\$150**
- Yellow----- 101C08 **\$150**
- Custom Colours----- 101C17 **\$150**

**③ SEAT WIDTH**

- 15" Seat Width ----- 101SW15 **No Charge**
- 16" Seat Width ----- 101SW16 **No Charge**
- 17" Seat Width ----- 101SW17 **No Charge**
- 18" Seat Width ----- 101SW18 **No Charge**
- 19" Seat Width ----- 101SW19 **No Charge**
- A  20" Seat Width ----- 101SW20 **No Charge**
- A  21" Seat Width (Weight Limit 350 lbs)----- 101SW21HD1 **\$475**
- A,B  22" Seat Width (Weight Limit 350 lbs)----- 101SW22HD1 **\$475**
- A,B  23" Seat Width (Weight Limit 350 lbs)----- 101SW23HD1 **\$685**
- 24" Seat Width (Weight Limit 350 lbs)----- 101SW24HD1 **\$685**

**④ SEAT DEPTH**

- 15" Seat Depth ----- 101SD15 **No Charge**
- 16" Seat Depth ----- 101SD16 **No Charge**
- 17" Seat Depth ----- 101SD17 **No Charge**
- 18" Seat Depth ----- 101SD18 **No Charge**
- 19" Seat Depth ----- 101SD19 **No Charge**
- 20" Seat Depth ----- 101SD20 **\$110**
- 21" Seat Depth ----- 101SD21 **\$225**
- 22" Seat Depth ----- 101SD22 **\$225**
- 23" Seat Depth ----- 101SD23 **\$225**
- T  24" Seat Depth ----- 101SD24 **\$225**
- EZ Tilt Option** ----- **101AL4 No Charge**

**Standard Features:**

Tilt Function (45 Degrees)  
Seat Width 18"  
Seat Depth 18"  
Canes and Stroller Bar  
T-Style Height Adjustable Full Length Gel  
Seat to Floor 18"  
Rear Wheels 24"  
Casters 8"  
Multi Adjustable Axel Plate  
70 Degree Front Rigging w/aluminum flip up footplate  
Heel Loops  
Push to Lock Wheel Locks  
Aluminum Seat Pan  
Prism Ideal Cushion  
Adjustable Orion Back  
Two way Padded, Contoured Headrest  
Anti-tippers  
Choice of 2 colours - Black or Silver



**FOOTNOTES:**

A. Heavy Duty Package consisting of gusseted frame

B. Quick release axles are used on 23"-24" wide chairs.


R. Add 9.5" to seat width of chair on finished product (handrim to handrim)

T. EZ Tilt option required on seat depths of 18" or higher


**ORION II SEAT TO FLOOR HEIGHT MATRIX**

FRONT STF	CASTER	FORK HOLE	STEM LENGTH
13.25	4"	1	X-SHORT
14.5	5"	1	SHORT
15.0	5"	2	SHORT
15.5	5"	3	SHORT
16.0	6"	3	SHORT
16.5	7"	3	SHORT
16.0	5"	3	LONG
16.5	6"	3	LONG
17.0	7"	3	LONG
17.0	8"	3	LONG
18.0	8"	4	LONG
* 19.0	8"	4	LONG
** 20.0	8"	4	LONG


\* 1" RISER BLOCK REQUIRED  
\*\* 2" RISER BLOCK REQUIRED



1  
CASTER FORK  
X-SHORT



1  
2  
3  
CASTER FORK  
SHORT




1  
2  
3  
4  
CASTER FORK  
LONG

**⑤ SEAT-TO-FLOOR HEIGHT (up to 350 lbs weight capacity) - Orion II**

- C  13.25" Seat to Floor (20" Wheels and 4" Casters) ----- 101STF01 \$225
- C  13.25" Seat to Floor (16" Wheels and 4" Casters) ----- 101STF97 No Charge
- 14.5" Seat to Floor (20" rear and 5" casters) ----- 101STF02 No Charge
- 15" Seat to Floor (20" Wheels and 5" Casters) ----- 101STF03 No Charge
- D  15.5" Seat to Floor (22" Wheels and 5" Casters) ----- 101STF25 No Charge
- D  15.5" Seat to Floor (12" Wheels and 6" Casters) ----- 101STF26 \$195
- 16" Seat to Floor (20" Wheels and 5" Casters) ----- 101STF27 No Charge
- 16" Seat to Floor (12" Wheels and 6" Casters) ----- 101STF28 \$195
- 16" Seat to Floor (24" Wheels and 6" Casters) ----- 101STF29 No Charge
- D  16" Seat to Floor (22" Wheels and 6" Casters) ----- 101STF06 No Charge
- 16 1/2" Seat to Floor (24" Wheels and 6" Casters) ----- 101STF30 No Charge
- 16 1/2" Seat to Floor (24" Wheels and 7" Casters) ----- 101STF31 No Charge
- 16 1/2" Seat to Floor (22" Wheels and 7" Casters) ----- 101STF08 No Charge
- 17" Seat to Floor (24" Wheels and 7" Casters) ----- 101STF32 No Charge
- 17" Seat to Floor (22" Wheels and 8" Casters) ----- 101STF33 No Charge
- 18" Seat to Floor (24" Wheels and 8" Casters) ----- 101STF11 Standard
- 19" Seat to Floor (24" Wheels and 8" Casters) 1" Spacer ----- 101STF34 No Charge
- 20" Seat to Floor (24" Wheels and 8" Casters) 2" Spacer ----- 101STF35 \$110

SUPER LOW

FRONT STF	REAR WHEEL	AXEL POS.
13.25	20"	SUPERLOW
13.25	16"	2
14.5	20"	1
14.5	16"	3
15.0	20"	2
15.5	22"	1
15.5	12"	4
16.0	24"	1
16.0	22"	2
16.0	20"	3
16.0	12"	4
16.5	24"	1
16.5	22"	2
16.5	24"	1
17.0	24"	2
17.0	22"	3
18.0	16"	4
18.0	24"	3
19.0	24"	3
20.0	24"	3



\* 1" RISER BLOCK REQUIRED  
\*\* 2" RISER BLOCK REQUIRED

**⑥ REAR WHEEL RIM**

- Mag Wheels ----- 101RW09 No Charge
- Spoke Wheels ----- 101RW10 \$75

**⑦ REAR WHEEL TIRE**

- Airless (Solid Rubber) ----- 101RT1 Standard
- Pneumatic ----- 101RT2 No Charge

**⑧ HANDRIM OPTIONS**

- Chrome ----- 101HR1 Standard
- Plastic Coated ----- 101HR2 \$160
- Omit Handrims ----- 101HR3 less net \$15

**FOOTNOTES:**

- C. Only available with pull to lock brakes.
- D. With 12" rear wheel you can only have push to lock brakes.

**9 ARM TYPE**

- T-Style Adjustable Height ----- 101AR04 Standard
- T-Style Adjustable Height Pin Release - Youth ----- 101AR05 \$280
- T-Style Adjustable Height Pin Release - Adult ----- 101AR06 \$300
- T-Style Fixed Height ----- 101AR02 \$250
- Flip Back Adjustable Height ----- 101AR10 \$225
- 2-point Fixed Height Heavy Duty Arm ----- 101AR13 \$250

**10 ARM PAD**

- Desk length Vinyl 10" (pair)----- 101AP1 No Charge
- Full length Vinyl 14" (pair)----- 101AP2 No Charge
- Contour Gel 17" (pair)----- 101AP3 Standard
- Contour Gel Wide 16" (pair)----- 101AP4 \$150

**11 PUSH OPTIONS**

- Canes & Stroller Bars Fixed Low ----- 101SP1 Standard
- Canes & Stroller Bars Fixed High ----- 101SP2 No Charge
- Canes & Stroller Bars Fold Down Low ----- 101SP3 \$250
- Canes & Stroller Bars Fold Down High ----- 101SP4 \$250
- Stroller Bars Fixed Low----- 101SH1 \$150
- Stroller Bars Fixed High----- 101SH2 \$150
- Stroller Bars Adjustable Angle Fold Down Low----- 101SH3 \$250
- Stroller Bars Adjustable Angle Fold Down High----- 101SH4 \$250
- Stroller Bars Adjustable Hinged (16"-20" Width Adj.)----- 101SH5 \$300
- Stroller Bars Adjustable Hinged (20"-24" Width Adj.)----- 101SH6 \$300
- Canes & Stroller Bars Recline Low----- 101SP5 \$500
- Canes & Stroller Bars Recline High----- 101SP6 \$500
- Stroller Bar Recline Low----- 101SP11 \$500
- Stroller Bar Recline High----- 101SP12 \$500
- Stroller Bar Fixed Dual Bend ----- 101SP15 No Charge
- Back Cane Angle (83, 90, 97, 104,111) Specify \_\_\_\_\_ 101BA1 No Charge

**12 TILT TRIGGER LOCATION**

- Trigger Location - Push Cane----- 101TH1 No Charge  
Side Trigger Handle not available with Adjustable Hinge Stroller.
- Trigger Location - Stroller Bar----- 101TH2 No Charge
- Foot Tilt Release (Mounted onto left side) ----- 101TH3 \$95
- Foot Tilt Release (Mounted onto right side) ----- 101TH4 \$95

**13 BACK OPTIONS (CHOOSE ONE)**

- Orion Back ----- 101BK00 No Charge
- Prism Basic Back ----- 101BK01 \$150
- Prism Ultra Back ----- 101BK06 \$250
- Prism Air Back ----- 101BK11 \$400
- Prism TrueFitt Back ----- 101BK16 \$650
- Prism Max Ultra Back ----- 101BK21 \$300
- Prism Max Air Back ----- 101BK26 \$450
- Prism Max TrueFitt Back ----- 101BK31 \$700
- Lunar Back Mid Basic Lateral ----- 101BK39 \$550

Specify Back Size Height: 16"\_\_ 18"\_\_ 20"\_\_

For other Backs, please fill in the required Backs Order Form and attach it to this order. [Click here to see all Backs](#)

**14 BACK PINS**

- Longer Pins for Back Hardware 1" Growth ----- 101BP1 \$50
- Longer Pins for Back Hardware 2" Growth ----- 101BP2 \$50

**FOOTNOTES:**

- E. Height adjustability of 4" and depth adjustability of 1"
- F. Max Backs come in one of the following colours: Black, Navy Blue, Royal Blue, Forest Green or Burgundy
- G. Please refer to Backs order form for for stock, non-stock, bariatric and custom sizes

**15 SEAT OPTIONS (CHOOSE ONE)**

<input type="checkbox"/> Prism Ideal Cushion -----	101SC01	No Charge
<input type="checkbox"/> Prism Elite Cushion -----	101SC04	\$250
<input type="checkbox"/> Prism Supreme Cushion - -----	101SC07	\$375
<input type="checkbox"/> Prism Supreme II Plus Cushion -----	101SC10	\$375
<input type="checkbox"/> Prism Supreme Air Cushion -----	101SC13	\$500
<input type="checkbox"/> Omit Cushion -----	101SC20	No Charge

Standard Size Cushion Depth: 16", 18" and 20"

For other Cushions, please fill in the required Cushions Order Form and attach it to this order. [Click here to see all Cushions](#)

**16 FRONT RIGGING**

<input type="checkbox"/> 70 Degree Pin Style -----	101FR01	Standard
<input type="checkbox"/> 80 Degree Pin Style -----	101FR02	\$135
<input type="checkbox"/> 60 Degree Swingaway Upgraded -----	101FR03	\$215
<input type="checkbox"/> 70 Degree Swingaway Upgraded -----	101FR04	\$215
<input type="checkbox"/> 90 Degree Swingaway Upgraded -----	101FR05	\$230
H <input type="checkbox"/> Elevating Legrest (pair) -----	101FR10	\$550
<input type="checkbox"/> Omit Front Rigging -----	101FR11	less net \$20
<input type="checkbox"/> Omit Footrest Hanger -----	101FR12	less net \$10

**17 FOOTPLATE TYPE**

<input type="checkbox"/> Composite Flip Up Footplates Small-----	101FP01	No Charge
<input type="checkbox"/> Composite Flip Up Footplates Large-----	101FP02	No Charge
<input type="checkbox"/> Aluminum Flip Up Footplates Small-----	101FP05	Standard
<input type="checkbox"/> Aluminum Flip Up Footplates Medium-----	101FP06	Standard
<input type="checkbox"/> Aluminum Flip Up Footplates Large-----	101FP07	Standard
<input type="checkbox"/> Composite Angle Adj. Flip-Up Footplates Medium (18" UP) ----	101FP09	\$120
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock Large-----	101FP16	\$150
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock X-Large-----	101FP17	\$170
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Small-----	101FP10	\$120
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Medium-----	101FP11	\$120
<input type="checkbox"/> Aluminum Angle Adj Flip Up Footplates Large-----	101FP12	\$120
<input type="checkbox"/> Composite Angle Adj. Flip-Up Footplates Spring Return Lg-----	101FP04	\$150
I <input type="checkbox"/> One Piece Flip-Up Footboard, Aluminum (15"-20") -----	101FP19	\$160
I <input type="checkbox"/> One Piece Flip-Up Footboard, Aluminum (21"-24" Custom)-----	101FP20	\$195
<input type="checkbox"/> One Piece Flip-Up Footboard, Padded -----	101FP26	\$210

**18 FOOTREST /LEGREST ACCESSORIES**

<input type="checkbox"/> Heel Loops -----	101FA01	Standard
<input type="checkbox"/> Omit Heel Loops -----	101FA02	less net \$7
<input type="checkbox"/> Padded Calf Strap -----	101FA03	\$70
J <input type="checkbox"/> Swingway Calf Pad (pair) -----	101FA04	\$185
<input type="checkbox"/> Foam Guard on Elevating Legrest (pair) -----	101FA05	\$70
<input type="checkbox"/> Foam Pad One Piece Footboard -----	101FA06	\$50
<input type="checkbox"/> Foam Pad on Footplate (pair) -----	101FA07	\$70




**19 WHEEL LOCKS**

<input type="checkbox"/> Push-to-Lock Wheel Lock -----	101WL1	Standard
<input type="checkbox"/> Pull-to-Lock Wheel Lock -----	101WL2	\$40
<input type="checkbox"/> Wheel Lock Extension -----	101WL3	No Charge
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) RT -----	101WL8	\$356
<input type="checkbox"/> Attendant Wheel Lock (Foot Activated) LT -----	101WL9	\$356
<input type="checkbox"/> Berner Push Wheel Lock -----	101WL14	\$95
<input type="checkbox"/> Berner Pull Wheel Lock -----	101WL15	\$95


**FOOTNOTES:**

- H. Includes Swingway Calf Pads
- I. Not available with 7" or 8" casters when footboard is facing back on 90 degree swing away

**COMPOSITE FOOTPLATES**




SMALL MEDIUM



**ALUMINUM FOOTPLATES**

ADULT ALUMINUM ANGLE ADJUSTABLE FOOTPLATE  
SMALL MEDIUM LARGE



**20 HEADREST HARDWARE OPTIONS**

<input type="checkbox"/>	Fixed Headrest Hardware -----	101HE01	No Charge
<input type="checkbox"/>	Fixed Headrest Curved Hardware -----	101HE15	No Charge
<input type="checkbox"/>	Fixed Headrest Hardware (Plus 4") -----	101HE02	\$20
<input type="checkbox"/>	Fixed Headrest Curved Hardware (Plus 4") -----	101HE16	\$20
<input type="checkbox"/>	Single Multi Axial Headrest Hardware -----	101HE03	\$158
<input type="checkbox"/>	Single Multi Axial Curved Headrest Hardware -----	101HE17	\$175
<input type="checkbox"/>	Single Multi Axial Headrest Hardware (+ 4") -----	101HE04	\$183
<input type="checkbox"/>	Single Multi Axial Curved Headrest Hardware (Plus 4") -----	101HE18	\$200
<input type="checkbox"/>	Dual Multi Axial Headrest Hardware -----	101HE05	\$195
<input type="checkbox"/>	Dual Multi Axial Curved Headrest Hardware -----	101HE19	\$225
<input type="checkbox"/>	Dual Multi Axial Headrest Hardware (Plus 4") -----	101HE06	\$220
<input type="checkbox"/>	Dual Multi Axial Curved Headrest Hardware (Plus 4") -----	101HE20	\$250
<input type="checkbox"/>	Infinite Multi Axial Headrest Hardware -----	101HE21	\$225
<input type="checkbox"/>	Infinite Multi Axial Curved Headrest Hardware -----	101HE22	\$250
<input type="checkbox"/>	Infinite Multi Axial Hardware (Plus 4") -----	101HE23	\$245
<input type="checkbox"/>	Infinite Multi Axial Curved Hardware (Plus 4") -----	101HE24	\$270
<input type="checkbox"/>	Links Multi Adjustable Headrest (Comes with 2 links) -----	101HE10	\$200
<input type="checkbox"/>	Links Multi Adjustable Curved Headrest (Comes with 2 links) --	101HE11	\$240
<input type="checkbox"/>	Additional Link for Headrest -----	101HE12	\$50
<input type="checkbox"/>	2 Additional Links for Headrest -----	101HE13	\$100
<input type="checkbox"/>	3 Additional Links for Headrest -----	101HE14	\$150

**21 HEADREST ADAPTOR BRACKETS**

<input type="checkbox"/>	2 - Hole Headrest Adaptor Bracket with Clamp -----	101HEBA2	No Charge
<input type="checkbox"/>	Universal Headrest Adaptor Bracket with Clamp -----	101HEBA	\$61

**22 HEADREST FOAM PAD OPTIONS**

<input type="checkbox"/>	Small Headrest Foam Pad -----	101HF1	No Charge
<input type="checkbox"/>	Small Headrest with Wings Foam Pad -----	101HF2	\$90
<input type="checkbox"/>	Medium Headrest Foam Pad -----	101HF7	No Charge
<input type="checkbox"/>	Medium Headrest with Wings Foam Pad -----	101HF8	\$90
<input type="checkbox"/>	Large Headrest Foam Pad -----	101HF3	No Charge
<input type="checkbox"/>	Large Headrest with Wings Foam Pad -----	101HF4	\$95
<input type="checkbox"/>	Pro-Fit Small Headrest -----	101HF12	\$100
<input type="checkbox"/>	Pro-Fit Medium Headrest -----	101HF13	\$105
<input type="checkbox"/>	Pro-Fit Large Headrest -----	101HF14	\$110
<input type="checkbox"/>	Omit Headrest -----	101HE09	No Charge

**23 ANTI TIP**

<input type="checkbox"/>	Quick Flip Anti-Tippers -----	101MA23	Standard
<input type="checkbox"/>	Fixed Anti-Tipper -----	101MA22	No Charge

**24 MISCELLANEOUS ACCESSORIES**

<input type="checkbox"/>	Tilt Locking Collar -----	101MA01	\$20
<input type="checkbox"/>	Recline Locking Collar -----	101MA02	\$20
<input type="checkbox"/>	Auto Seat Positioning Strap -----	101MA03	\$52
<input type="checkbox"/>	Airline Seat Positioning Strap -----	101MA04	\$76
<input type="checkbox"/>	Velcro Release Positioning Strap -----	101MA05	\$48
<input type="checkbox"/>	Foam Pad on Arm Guard (pair) -----	101MA06	\$140
<input type="checkbox"/>	Padded Lap Tray 15"-16" -----	101MA09	\$285
<input type="checkbox"/>	Padded Lap Tray 17"-18" -----	101MA10	\$285
<input type="checkbox"/>	Padded Lap Tray 19"-20" -----	101MA11	\$285
<input type="checkbox"/>	Padded Lap Tray 21"-22" -----	101MA12	\$385
<input type="checkbox"/>	Padded Lap Tray 23"-24" -----	101MA13	\$385
<input type="checkbox"/>	Half Padded Lap Tray, left side -----	101MA14	\$160
<input type="checkbox"/>	Half Padded Lap Tray, right side -----	101MA15	\$160
<input type="checkbox"/>	Elbow Stop, left -----	101MA16	\$94
<input type="checkbox"/>	Elbow Stop, right -----	101MA17	\$94
<input type="checkbox"/>	Shoulder Lateral Pad (each) -----	101MA18	\$90
<input type="checkbox"/>	Adjustable Side Pad -----	101MA19	\$157
<input type="checkbox"/>	IV Pole, telescopic -----	101MA20	\$168
<input type="checkbox"/>	Oxygen Tank Holder -----	101MA21	\$173
<input type="checkbox"/>	Quick Release Axle -----	101AL1	\$110
<input type="checkbox"/>	Adductor Pad 10W x 4.5H w/Flush Hooks (pair) -----	101MA24	\$175
<input type="checkbox"/>	Adductor Pad 10W x 4.5H w/Offset Hooks (pair) -----	101MA25	\$175
<input type="checkbox"/>	Adductor Pad 12W x 4.5H w/Flush Hooks (pair) -----	101MA26	\$175
<input type="checkbox"/>	Adductor Pad 12W x 4.5H w/Offset Hooks (pair) -----	101MA27	\$175
<input type="checkbox"/>	Adductor Pad 14W x 4.5H w/Flush Hooks (pair) -----	101MA28	\$175
<input type="checkbox"/>	Adductor Pad 14W x 4.5H w/Offset Hooks (pair) -----	101MA29	\$175
<input type="checkbox"/>	Spoke Guards available on 20", 22" and 24" wheels -----	101MA49	\$130
<input type="checkbox"/>	Orion II Front and Rear Transit Hooks -----	101MA30	\$250

K  
L  
N

**FOOTNOTES:**

- K. Only available with Stroller Bar push option.
- L. Not available with Flip Back Arms
- M. Weight Capacity of Orion II with Power Tilt is 250 lbs. Power Tilt can only be placed on Orion II chairs with seat width of 15" to 20" Power Tilt provides approximately 42 degrees of infinite tilt (3° to 45°) Includes: 2 Button Hand Control,
- N. Transit Hooks can only be used on Orion II wheelchairs up to 350 lbs weight capacity. Tested in accordance with ANSI/RESNA WC/19

<b>WARRANTY</b>
Future Mobility Healthcare warrants the Orion II to be free from defects in materials and workmanship for a period of two (2) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user. The batteries required for Power Tilt are warranted for six (6) months while the control box, battery charger, actuator and hand control are warranted for one (1) year.