


ADJUSTABLE TENSION BACK



The Adjustable Tension Back™ was designed to provide additional comfort and adjustability. With two fixed straps and two removable velcro straps, the user is able to apply the tension where it is needed most. This allows for all day comfort and pressure relief.

Prism Series Backs STYLE: Adjustable Tension Back				Tension Back High - Shipment within 3 business days			
Tension Back Low - Same Day Shipment				Part Number	Size W x H	Part Number	Size W x H
ADP: WAMD - Backs under 18" Wide				TB1416	14" x 16"	TB1816	18" x 16"
WAME - Backs over 18" Wide				TB1516	15" x 16"	TB1916	19" x 16"
Part Number	Size W x H	Part Number	Size W x H	TB1616	16" x 16"	TB2016	20" x 16"
TB1413	14" x 13"	TB1813	18" x 13"	TB1716	17" x 16"		
TB1513	15" x 13"	TB1913	19" x 13"	Custom Sizes - Shipment within 7 business days			
TB1613	16" x 13"	TB2013	20" x 13"	Part Number	Size W x H		
TB1713	17" x 13"			TBCUS	Custom		



Product Highlights:

- ✓ Attractive design
- ✓ Pocket for personal items
- ✓ Two fixed straps and two removable straps
- ✓ All day comfort and pressure relief
- ✓ **Product Weight:** 2.6 lbs (18"W x 18"H)
- User Weight:** 250 lbs
- Warranty:** 1 year limited