



**FUTURE**  
MOBILITY HEALTHCARE INC.



BK PHANTOM PB 2830 REV02

# **PARTS MANUAL**

FUTURE MOBILITY HEALTHCARE APPRECIATES YOUR BUSINESS.  
OUR CUSTOMER SERVICE REPRESENTATIVES CAN BE CONTACTED AT:

**TO PLACE ORDERS:  
(905) 671-1661**

FOR YOUR CONVENIENCE YOU MAY ALSO PLACE ORDERS BY FAX AT:

**FAX: (905) 671-3377**

PLEASE VISIT OUR WEBSITE AT:

**[www.futuremobility.ca](http://www.futuremobility.ca)**

\*Part prices and configurations included in this Parts Manual are subject to change without prior notice.

\* All prices are shown as per unit individual price.


## Contents

PHANTOM BACK .....	5
<b>PHANTOM BACK HARDWARE .....</b>	<b>5</b>



# PHANTOM BACK

## PHANTOM BACK HARDWARE



## PHATON BACK HARDWARE TABLE

ITEM	NOTES	QTY	PART NO.	DESCRIPTION	PRICE/ UNIT
1		4	1000401	WASHER FLAT 1/4", ID .257, OD 0.50, THK 0.062	\$1.50
2		1	REF	GENERAL PHB SHELL	-
3		2	1002903	PH CLAMP MOUNT ST	\$45.00
4		2	1003156	PH POST CLAMP	\$75.00
5		2	1004037	SLEEVE BUSHING, PHANTOM BACK CLAMP -3/4"	\$20.00
6		2	1004068	2-317 N70 O RING 1- 5/16 OD	\$5.00
7		6	1004069	SCREW,SOC BUTTON HEAD 1/4-20 X 0.375	\$3.00
8		2	1004070	SCREW ,BUTTON HEAD 1/4-20 X 5/8	\$3.00
		1	<b>1004923</b>	<b>PHANTOM BACK HARDWARE KIT (ITEMS 1, 3-8)</b>	<b>\$170.00</b>

NOTES: