

ORDER FORM

Magellan



TAG:

ACCOUNT INFORMATION

Acc. # :

Contact Name:

Company Name:

SHIPPING INFORMATION

PO #:

Street:

City:

Province:

Postal Code:

Country:



Please call Customer Service if you require any unique customization or if you have inquiries
1-888-737-4011

1 FRAME TYPE

201 **MAGELLAN - WPABFM005** \$6045 X

CONTROLLER CUSTOMIZATION :

MOUNTING SIDE: R L	MOUNTING OPTIONS: SWING-AWAY	FIXED
201CS1 201CS2	201JM2	201JM1
JOYSTICK OPTIONS: VR2 PM90 +\$0	RNET PM80 +\$1345	STRAIGHT EXTENSION
201JY1	201JY4	201JY2 +\$58
		T EXTENSION 201JY3 +\$83

2 FRAME FINISH

201C01	MAT BLACK	_____	\$0 X
201C02	SILVER LINING	_____	\$0 X
201C03	SKY BLUE		\$0X
201C04	HARLEQUIN GREEN		\$0X
201C05	PUMPKIN ORANGE		\$0X
201C06	AMETHYST PURPLE		\$0X
201C07	IMPERIAL RED		\$0X
201C08	BUMBLEBEE YELLOW		\$0X
201C07	SILVER VAIN		\$0X
201C08	BUBBLEGUM PINK		\$0X

Standard Features:

- SEAT WIDTH - 16" - 24"
- SEAT WIDTH INCLUDING JOYSTICK - 22 1/4" - 30 1/4"
- BASE WIDTH - 24 1/4"
- LENGTH W. FOOT BOARD IN USE - 37"
- LENGTH W. SWING-AWAY FOOT RESTS IN USE - 40" - 45"
- HEIGHT - 40" - 45"
- SEAT TO FLOOR - 17" - 19"
- REAR WHEEL - 12", 16", 20", 22" OR 24" COMPOSITE URETHANE
- SEAT DEPTH - 15" - 24"
- OVERALL WIDTH - NARROW SEAT WIDTH + 9.5", REGULAR SEAT WIDTH + 10.0"
- FRONT & REAR CASTERS - 6" X 2 SEMI-PNEUMATIC
- FRONT RIGGING - SWINGAWAY, DUAL SWINGAWAY OR DUAL ELEVATING LEG RESTS
- ADJUSTABLE BACK ANGLE - -7, 0, 7, 14, 21 DEGREES
- BACK STYLE - LOW PUSH CANES WITH STROLLER BAR. HIGH PUSH CANES WITH STROLLER BAR. 16", 18", 20" ORION BACK
- GROUND CLEARANCE - 3"
- TURNING RADIUS - 20.3"
- RANGE PER CHARGE - 13.4 MILES (21.6KM)
- MAX SPEED - 4.3MPH (7KM/H)
- ARM STYLE - ADJUSTABLE HEIGHT T-STYLE, ADJUSTABLE HEIGHT PIN RELEASE T-STYLE, ADJUSTABLE HEIGHT FLIP BACK
- PRODUCT WEIGHT - 200 LBS (WITH BATTERIES) OR 132LBS (WITHOUT BATTERIES)
- WEIGHT CAPACITY - 300 LBS (136 KG)

3 DRIVE WHEELS					
201DT2	AIRLESS	\$0 X	201DT1	PNEUMATIC	\$145 X
4 CHARGERS					
201CH2	6A ON BOARD	\$0 X	201CH3	8A OFF BOARD	\$329 X
5 BATTERIES					
	201BT2	8G24 BATTERIES (2 PER CHAIR)		\$500 X	
6 STANDARD SEATING 201SE6					
		\$842			
		NO	YES		
(INCLUDES SEAT FRAME, SOLID SEAT PAN, 2-POINT ADJUSTABLE HEIGHT FLIPBACK REMOVABLE ARMS)					
7 SEAT WIDTH					
201SW16	16" WIDTH	\$0 X	201SW21	21" WIDTH	\$0 X
201SW17	17" WIDTH	\$0 X	201SW19	22" WIDTH	\$274 X
201SW18	18" WIDTH	\$0 X	201SW20	23" WIDTH	\$274 X
201SW19	19" WIDTH	\$274 X	201SW20	24" WIDTH	\$274 X
201SW20	20" WIDTH	\$274 X			
8 SEAT DEPTH					
201SD17	17" SEAT DEPTH	\$0 X	201SD19	19" SEAT	\$0 X
201SD18	18" SEAT DEPTH	\$0 X	201SD20	20" DEPTH	\$0 X
9 CANE OPTIONS					
201PC3	PUSH CANE FIXED	\$0 X	201SH5	PUSH CANE MANUAL RECLINE	\$0 X
10 BACK OPTIONS (CONTINUED ON NEXT PAGE)					
201BK48	ORION BACK (16", 18", 20") AVAIL- ABLE HEIGHT ONLY	\$0 X	201BK11	PRISM AIR BACK ON CHAIR STOCK	\$400 X
201BK01	PRISM BASIC BACK ON CHAIR STOCK	\$150 X	201BK12	PRISM AIR BACK ON CHAIR NON STOCK	\$475 X
201BK02	PRISM BASIC BACK ON CHAIR NON STOCK	\$175 X	201BK13	PRISM AIR BACK ON CHAIR BARIAT- RIC	\$600 X
201BK03	PRISM BASIC BACK ON CHAIR BARIATRIC	\$250 X	201BK14	PRISM AIR BACK CUSTOM ON CHAIR	\$1075 X
201BK04	PRISM BASIC BACK ON CHAIR CUSTOM	\$825 X	201BK15	PRISM AIR BACK BARIATRIC CUSTOM	\$1375 X
201BK05	PRISM BASIC BACK BARIATRIC CUSTOM ON CHAIR	\$1,075 X	201BK16	PRISM TRUEFIT BACK ON CHAIR STOCK	\$650 X
201BK06	PRISM ULTRA BACK ON CHAIR STOCK	\$250 X	201BK17	PRISM TRUEFIT BACK ON CHAIR NON STOCK	\$700 X
201BK07	PRISM ULTRA BACK ON CHAIR NON STOCK	\$300 X	201BK18	PRISM TRUEFIT BACK ON CHAIR	\$950 X
201BK08	PRISM ULTRA BACK ON CHAIR BARIATRIC	\$350 X	201BK19	PRISM TRUEFIT BACK CUSTOM ON CHAIR	\$1025 X
201BK09	PRISM ULTRA BACK CUSTOM ON CHAIR	\$1010 X	201BK20	PRISM TRUEFIT BACK BARIATRIC CUSTOM ON CHAIR	\$1350 X
201BK10	PRISM ULTRA BACK BARIATRIC CUSTOM ON CHAIR	\$1260 X	201BK21	MAX ULTRA BACK ON CHAIR STOCK	\$300 X

BACK OPTIONS (CONTINUED)

201BK22	MAX ULTRA BACK ON CHAIR NON STOCK	\$350 X	201BK34	MAX TRUFITT BACK CUSTOM ON CHAIR	\$1160 X
201BK23	MAX ULTRA BACK ON CHAIR BARIATRIC	\$400 X	201BK35	MAX TRUFITT BACK BARIATRIC CUSTOM ON CHAIR	\$1435 X
201BK24	MAX ULTRA BACK CUSTOM ON CHAIR	\$1095 X	201BK36	LUNAR BACK LOW NO LATERAL ON CHAIR	\$400 X
201BK25	MAX ULTRA BACK BARIATRIC CUSTOM ON CHAIR	\$1345 X	201BK37	LUNAR BACK CUSTOM ON CHAIR	\$475 X
201BK26	MAX AIR BACK ON CHAIR STOCK	\$450X	201BK38	LUNAR BACK MID NO LATERAL ON CHAIR	\$490 X
201BK27	MAX AIR BACK ON CHAIR NON STOCK	\$525 X	201BK39	LUNAR BACK MID BASIC LATERAL ON CHAIR	\$550 X
201BK28	MAX AIR BACK ON CHAIR BARIATRIC	\$650 X	201BK40	LUNAR BACK MID ULTRA LATERAL ON CHAIR	\$610 X
201BK29	MAX AIR BACK CUSTOM ON CHAIR	\$1160 X	201BK41	LUNAR BACK FULL NO LATERAL ON CHAIR	\$510 X
201BK30	MAX AIR BACK BARIATRIC CUSTOM ON CHAIR	\$1410 X	201BK42	LUNAR BACK FULL BASIC LATERAL ON CHAIR	\$570 X
201BK31	MAX TRUFITT BACK ON CHAIR STOCK	\$700 X	201BK43	LUNAR BACK FULL ULTRA LATERAL ON CHAIR	\$630 X
201BK32	MAX TRUFITT BACK ON CHAIR NON STOCK	\$750 X	201BK44	LUNAR BACK CUSTOM ON CHAIR	\$950 X
201BK33	MAX TRUFITT BACK ON CHAIR BARIATRIC	\$1000 X			

201BK49 OMIT BACK - \$30 X

SPECIFY BACK SIZE

BACK PINS

201BP1	LONGER PINS FOR BACK HARDWARE 1" GROWTH	\$50 X	201BP2	LONGER PINS FOR BACK HARDWARE 2" GROWTH	\$50 X
--------	-----------------------------------------	--------	--------	-----------------------------------------	--------

CUSHION OPTIONS

201SE1	SEAT PAN	\$0 X	201SC08	PRISM SUPREME AIR CUSHION - CUSTOM SIZES	\$400 X
201SC01	PRISM IDEAL CUSHION - STOCK SIZES	\$0 X	201SC09	PRISM SUPREME CUSHION - CUSTOM SIZES	\$685 X
201SC02	PRISM IDEAL CUSHION - NON-STOCK SIZES	\$0 X	201SC10	PRISM SUPREME II PLUS CUSHION STOCK	\$375 X
201SC03	PRISM IDEAL CUSHION - CUSTOM SIZES	\$0 X	201SC11	PRISM SUPREME II PLUS CUSHION NON-STOCK	\$400 X
201SC04	PRISM ELITE CUSHION - NON-STOCK SIZES	\$250 X	201SC12	PRISM SUPREME II PLUS CUSHION CUSTOM	\$685 X
201SC05	PRISM ELITE CUSHION - CUSTOM SIZES	\$275 X	201SC13	PRISM SUPREME AIR CUSHION - STOCK SIZES	\$500 X
201SC06	PRISM SUPREME CUSHION - STOCK SIZES	\$425 X	201SC14	PRISM SUPREME AIR CUSHION - NON-STOCK SIZES	\$575 X
201SC07	PRISM SUPREME CUSHION - NON-STOCK SIZES	\$375 X	201SC15	PRISM SUPREME AIR CUSHION CUSTOM ON CHAIR	\$785X

201BK49 OMIT CUSHION - \$10 X

ARM PADS

201AR24	OVAL FULL LENGTH, RIGHT	\$0 X	201AR27	WIDE MULTIPOSITION FULL LENGTH, LEFT	\$98 X
201AR25	OVAL FULL LENGTH, LEFT	\$0 X	201AP2	CONTOUR GEL WIDE 16" (PAIR)	\$150 X
201AR26	WIDE MULTIPOSITION FULL LENGTH, RIGHT	\$98 X			

BACK ANGLE

201BA2	FIXED (+3 TO -20 DEGREES)	\$0 X	201BA3	MANUAL RECLINE (+5 TO -40 DEGREES)	\$365 X
--------	---------------------------	-------	--------	------------------------------------	---------

FRONT RIGGING

201FR06	60 DEGREE SWINGAWAY UPGRAD-ED	\$150 X	201FR01	70 DEGREE PIN STYLE	\$135 X
201FR07	70 DEGREE SWINGAWAY UPGRAD-ED	\$150 X	201FR02	80 DEGREE PIN STYLE	\$135 X
201FR09	90 DEGREE SWINGAWAY UPGRAD-ED	\$175 X	201FR11	OMIT FRONT RIGGING	-\$10 X
201FR10	UPGRADED ELEVATING LEGREST	\$361 X	201FP21	ONE PIECE FOOTBOARD	\$0 X

FOOT PLATE TYPE

201FP09	COMPOSITE ANGLE ADJ. FLIP-UP FOOT-PLATES	\$113 X	201FP04	COMPOSITE ADJ ANGLE FLIP-UP SPRING RETURN FOOTPLATES	\$150 X
201FP11	ALUMINUM ANGLE ADJ FLIP UP FOOT-PLATES	\$113 X	201FP16	COMPOSITE ANGLE ADJ. FLIP-UP CONE LOCK LARGE	\$150 X
201FP02	COMPOSITE-UP FOOTPLATES	\$0 X	201FP17	COMPOSITE ANGLE ADJ. FLIP-UP CONE LOCK X-LARGE	\$170 X

FOOT RESTS & ACCESSORIES

201FA01	HEEL LOOPS	\$0 X	201FA06	FOAM PAD ONE PIECE FOOTBOARD	\$50 X
201FA02	OMIT HEEL LOOPS	-\$7 X	201FA07	FOAM PAD ON FOOTPLATE (PAIR)	\$70 X
201FA05	FOAM GUARD ON ELEVATING LEGREST (PAIR)	\$70 X			

HEAD REST HARDWARE OPTIONS

201HE01	FIXED	\$0 X	201HE11	LINKS MULTI ADJUSTABLE CURVED HEADREST	\$183 X
201HE03	SINGLE MULTI AXIAL CURVED	\$0 X	201HE12	ADDITIONAL LINK FOR HEADREST	\$200 X
201HE05	DUAL MULTI-AXIAL HEADREST	\$20 X	201HE13	2 ADDITIONAL LINKS FOR HEAD-REST	\$195 X
201HE02	FIXED HEADREST HARDWARE ASIST* (PLUS 4")	\$20 X	201HE14	3 ADDITIONAL LINKS FOR HEAD-REST	\$225 X
201HE04	SINGLE MULTI-AXIAL HEADREST HARDWARE (PLUS 4")	\$158 X	201HEBA	UHRBA	\$220 X
201HE06	DUAL MULTI-AXIAL HEADREST HARDWARE (PLUS 4")	\$175 X			
201HE10	LINKS MULTI ADJUSTABLE HEAD-REST	\$225 X			

HEAD REST FOAM PADS


201HF1	SMALL FOAM PAD	\$0 X	201HF4	LARGE W. WINGS	\$95 X
201HF2	SMALL PAD W. WINGS	\$90 X	201HF8	MEDIUM PAD W. WINGS	\$90 X
201HF7	MEDIUM FOAM PAD	\$0 X	201HF3	LARGE FOAM PAD	\$0 X
201HE09	OMIT HEAD REST	\$0 X			

MISC. ACCESSORIES

201MA32	CHEST BELT	\$106 X	201MA57	OMIT LAP BELT	\$0 X
201MA30	ANCHOR BRACKETS	\$60 X			

Warranty:

Base frame, seat frame and wheel hubs (not tires) 60 months limited warranty. Controller, charger and cables 24 months limited warranty. Motors and gearboxes 18 months limited warranty. Upholstery, plastic and rubber components, bearings, front rigging, forks and all other parts 6 months limited warranty.



**FUTURE MOBILITY
HEALTHCARE**