



Account Information		Order Date: _____	
Account #:	_____	P.O.#	_____
Business Name:	_____	Phone:	_____
Ship to Address:	_____	City:	_____
Fax:	_____	Province:	_____
Contact:	_____	Postal Code:	_____
TAG FOR:			

① FRAME TYPE

- Orion II 500- WMA5FM010 ----- 105 **\$2,700**

② FRAME FINISH

- Black ----- 105C01 No Charge
- Silver ----- 105C02 No Charge
- Custom Colours**
- Blue----- 105C03 \$150
- Green----- 105C04 \$150
- Orange----- 105C05 \$150
- Purple----- 105C06 \$150
- Red----- 105C07 \$150
- Yellow----- 105C08 \$150
- Custom Colour----- 105C17 \$150

Standard Features:

- Tilt Function (25 Degrees)
- Seat Width 24"-32"
- Seat Depth 15"-24"
- Fixed Stroller Bar
- 2-Point Heavy Fixed Arm
- Low Fixed Back
- Seat to Floor 16"
- Rear Wheels 20"
- Casters 6"
- Multi Adjustable Axle Plate
- 70 Degree Front Rigging w/footplate
- Heel Loops
- Spring Back Push to Lock Wheel Locks
- Aluminum Seat Pan
- Heavy Duty Back Upholstery
- Two way Padded, Contoured Headrest
- Anti-tippers
- Choice of 2 colours - Black or Silver

③ SEAT WIDTH

- ^A 24" Seat Width (Heavy Duty up to 500 lbs)----- 105SW24 \$1,200
- 25" Seat Width (Heavy Duty up to 500 lbs)----- 105SW25 \$1,200
- 26" Seat Width (Heavy Duty up to 500 lbs)----- 105SW26 \$1,200
- 27" Seat Width (Heavy Duty up to 500 lbs)----- 105SW27 \$1,200
- 28" Seat Width (Heavy Duty up to 500 lbs)----- 105SW28 \$1,200
- 29" Seat Width (Heavy Duty up to 500 lbs)----- 105SW29 \$1,200
- 30" Seat Width (Heavy Duty up to 500 lbs)----- 105SW30 \$1,200
- 31" Seat Width (Heavy Duty up to 500 lbs)----- 105SW31 \$1,200
- 32" Seat Width (Heavy Duty up to 500 lbs)----- 105SW32 \$1,200

④ SEAT DEPTH


- 15" Seat Depth ----- 105SD15 No Charge
- 16" Seat Depth ----- 105SD16 No Charge
- 17" Seat Depth ----- 105SD17 No Charge
- 18" Seat Depth ----- 105SD18 No Charge
- 19" Seat Depth ----- 105SD19 No Charge
- 20" Seat Depth ----- 105SD20 \$110
- 21" Seat Depth ----- 105SD21 \$225
- 22" Seat Depth ----- 105SD22 \$225
- 23" Seat Depth ----- 105SD23 \$225
- 24" Seat Depth ----- 105SD24 \$225



FOOTNOTES:

A. Double gas cylinder used on widths 24-32". Tilt 25 degrees

ORION II 500 FRONT STF HEIGHT AVAILABILITY									
REAR WHEEL		5" Castor							
12"					16.5"				
20"	15.5"								
22"	15.5"								
24"		16"	16.25"						
REAR WHEEL		6" Castor							
12"					16.5"				
20"		16"							
22"					16.5"				
24"					16.75"	17"	17.25"		
REAR WHEEL		7" Castor							
12"									
20"									
22"					16.5"	17"			
24"						17"	17.5"		
REAR WHEEL		8" Castor							
22"								17.75"	
24"								17.75"	18"



STF	REAR	CASTER	FORK	FORK	AXLE PLATE	STF
FRONT (in)	WHEEL (in)	SIZE (in)	SIZE	HOLE	HOLE	REAR (in)
16.5"	12"	6"	SHORT	C	D	15.25"
16.5"	12"	5"	SHORT	D	D	15.25"
15.5"	20"	5"	SHORT	B	B	14"
16"	20"	6"	SHORT	B	C	15"
15.5"	22"	5"	SHORT	B	A	14"
16.5"	22"	6"	LONG	B	B	15"
16.5"	22"	7"	LONG	B	B	15"
17"	22"	7"	SHORT	C	C	15.75"
17.75"	22"	8"	LONG	C	C	16"
16"	24"	5"	SHORT	C	A	14.75"
16.25"	24"	5"	LONG	C	A	14.75"
16.75"	24"	6"	LONG	C	A	14.75"
17.25"	24"	6"	LONG	D	B	15.75"
17"	24"	6"	SHORT	D	B	15.75"
17.5"	24"	7"	SHORT	D	B	15.75"
17"	24"	7"	SHORT	C	B	15.75"
17.75"	24"	8"	LONG	C	C	16.75"
18"	24"	8"	SHORT	D	C	16.75"

⑤ SEAT-TO-FLOOR HEIGHT

- 15.5" Seat to Floor (22" Wheels and 5" Casters) 105STF25 No Charge
- 15.5" Seat to Floor (20" Wheels and 5" Casters) 105STF04 No Charge
- 16" Seat to Floor (20" Wheels and 6" Casters) 105STF36 No Charge
- 16" Seat to Floor (24" Wheels and 5" Casters) 105STF37 No Charge
- 16.25" Seat to Floor (24" Wheels and 5" Casters) 105STF38 No Charge
- 16.5" Seat to Floor (12" Wheels and 5" Casters) 105STF39 \$195
- 16.5" Seat to Floor (12" Wheels and 6" Casters) 105STF40 \$195
- 16.5" Seat to Floor (22" Wheels and 6" Casters) 105STF07 No Charge
- 16.75" Seat to Floor (24" Wheels and 6" Casters) 105STF41 No Charge
- 17" Seat to Floor (22" Wheels and 7" Casters) 105STF09 No Charge
- 17" Seat to Floor (24" Wheels and 6" Casters) 105STF17 No Charge
- 17" Seat to Floor (24" Wheels and 7" Casters) 105STF32 No Charge
- 17.25" Seat to Floor (24" Wheels and 6" Casters) 105STF42 No Charge
- 17.5" Seat to Floor (24" Wheels and 7" Casters) 105STF19 No Charge
- 17.75" Seat to Floor (22" Wheels and 8" Casters) 105STF43 No Charge
- 17.75" Seat to Floor (24" Wheels and 8" Casters) 105STF44 No Charge
- 18.0" Seat to Floor (24" Wheels and 8" Casters) 105STF11 No Charge

FOOTNOTES:

B. With 12" rear wheel you can only use push to lock brakes.

⑥ REAR WHEEL RIM

- Mag Wheels 105RW09 No Charge
- Spoke Wheels 105RW10 \$75

⑦ REAR WHEEL TIRE

- Airless (Solid Rubber) ----- 105RT1 Standard
- Pneumatic ----- 105RT2 No Charge

⑧ HANDRIM OPTIONS

- Chrome ----- 105HR1 Standard
- Plastic Coated ----- 105HR2 \$160
- Omit Handrims----- 105HR3 less net \$15

⑨ ARM TYPE


- ^K 2-point Fixed Height Heavy Duty Arm ----- 105AR13 Standard
- T-Style Adjustable Height ----- 105AR04 \$280
- ^{K, O} Flip Back Fixed Height Heavy Duty Arm ----- 105AR14 \$350

FOOTNOTES:
K. Armpads are placed on inside part of chair to reduce interference with rear wheels
O. only available with vinyl arm pad

Armrest Matrix for T-Style Pin Release Arm						
Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range	Arm Pad Length L		
				Desk	Full	Curved
Adult	6.75"	7.50"	9.50" - 14.50"	10.25"	14.0"	17.0"
Youth	4.50"	5.0"	8.75" - 13.75"			

⑩ ARM PAD

- Desk length Vinyl 10" (pair)----- 105AP2 No Charge
- Full length Vinyl 14" (pair)----- 105AP1 No Charge
- Contour Gel 17" (pair)----- 105AP3 Standard
- Contour Gel Wide 16" (pair)----- 105AP4 \$150



⑪ PUSHBAR OPTIONS

- Stroller Bars Fixed Low ----- 105SP1 Standard
- Stroller Bars Fixed High ----- 105SP2 No Charge

⑫ TRIGGER LOCATION

- Trigger Location - Stroller Bar----- 105TH2 No Charge

⑬ BACK OPTIONS

- Heavy Duty Back Upholstery ----- 105BK50 Standard
- Orion Heavy Duty Back (16",18",20") available height ----- 105BK51 \$250
- Prism Basic Back Heavy Duty on chair ----- 105BK52 \$400
- Prism Basic Back Custom Heavy Duty on chair----- 105BK53 \$1,075
- Prism Ultra Back Heavy Duty on chair ----- 105BK54 \$500
- Prism Ultra Back Custom Heavy Duty on chair----- 105BK55 \$1,260
- Prism Air Back Heavy Duty on chair ----- 105BK56 \$600 *
- Prism Air Back Custom Heavy Duty on chair ----- 105BK57 \$1,325 *
- Prism TrueFitt Back Heavy Duty on chair ----- 105BK58 \$1,000
- Prism TrueFitt Back Custom Heavy Duty on chair----- 105BK59 \$1,350
- Prism Max Ultra Back Heavy Duty on chair ----- 105BK60 \$550
- Prism Max Ultra Back Custom Heavy Duty on chair ----- 105BK61 \$1,345
- Prism Max Air Back Heavy Duty on chair ----- 105BK62 \$800
- Prism Max Air Back Custom Heavy Duty on chair ----- 105BK63 \$1,410
- Prism Max TrueFitt Back Heavy Duty on chair ----- 105BK64 \$1,050
- Prism Max TrueFitt Back Custom Heavy Duty----- 105BK65 \$1,435
- ^M Lunar Back Heavy Duty Mid Basic Lateral on chair (22-26W)--- 105BK66 \$750
- ^M Lunar Back Heavy Duty Mid Basic Lateral on chair (28-32W)--- 105BK67 \$1,300
- ^M Lunar Back Heavy Duty Mid Ultra Lateral on chair (22-26W)----- 105BK68 \$810
- ^M Lunar Back Heavy Duty Mid Ultra Lateral on chair (22-26W)----- 105BK69 \$1,360
- ^N Lunar Back Heavy Duty Full Basic Lateral on chair (22-26W)----- 105BK70 \$800
- ^N Lunar Back Heavy Duty Full Basic Lateral on chair (28-32W)----- 105BK71 \$1,350
- ^N Lunar Back Heavy Duty Full Ultra Lateral on chair (22-26W)----- 105BK72 \$860
- ^N Lunar Back Heavy Duty Full Ultra Lateral on chair (28-32W)----- 105BK73 \$1,410
- Lunar Back Heavy Duty Custom on Chair (22W - 32W)----- 105BK74 \$1,800
- Omit Back ----- 105BK49 less net \$30
- ^D Specify Max Back Color: _____
- ^E Specify Back Size W____x H____

FOOTNOTES:
D. Max Backs come in one of the following colours: Black, Navy Blue, Royal Blue, Forest Green or Burgundy
E. Please refer to Backs order form for for stock, non-stock, heavy duty and custom sizes
L. Customer to specify trigger handle location for Orion II 500 wheelchairs only
M. Heights available 16H, 18H
N. Heights available 20H, 22H, 24H (refer to heavy duty lunar back order form for specific sizes)

14 BACK PINS (Prism Backs only)

<input type="checkbox"/> Longer Pins for Back Hardware 1" Growth-----	105BP2	\$50
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth-----	105BP3	\$50

15 SEAT OPTIONS

<input type="checkbox"/> Seat Pan -----	105SE1	Standard
<input type="checkbox"/> Prism Supreme Heavy Duty Cushion Stock on chair -----	105SC39	\$600
<input type="checkbox"/> Prism Supreme Heavy Duty Cushion Non Stock on chair -----	105SC40	\$625
<input type="checkbox"/> Prism Supreme Heavy Duty Cushion Custom -----	105SC41	\$985
<input type="checkbox"/> Omit Cushion -----	105SC20	less net \$10

16 FRONT RIGGING

<input type="checkbox"/> 70 Degree Pin Style -----	105FR01	Standard
<input type="checkbox"/> 80 Degree Pin Style -----	105FR02	\$135
<input type="checkbox"/> 60 Degree Swingaway Upgraded -----	105FR03	\$215
<input type="checkbox"/> 70 Degree Swingaway Upgraded -----	105FR04	\$215
F <input type="checkbox"/> Elevating Legrest (pair) -----	105FR10	\$550
<input type="checkbox"/> Omit Front Rigging -----	105FR11	less net \$20
<input type="checkbox"/> Omit Footrest Hanger -----	105FR12	less net \$10

18 FOOTPLATE TYPE

<input type="checkbox"/> Composite Flip-Up Footplates Large-----	105FP02	No Charge
<input type="checkbox"/> Aluminum Flip-Up Footplates Large-----	105FP07	Standard
<input type="checkbox"/> Aluminum Adjustable Angle Flip-Up Footplates Large-----	105FP12	\$120
<input type="checkbox"/> Composite Adjustable Angle Flip-Up Footplates Medium (18" up)	105FP09	\$120
<input type="checkbox"/> Composite Adjustable Angle Flip-Up Footplates Spring Return	105FP04	\$150
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock Large-----	105FP16	\$150
<input type="checkbox"/> Composite Angle Adj. Flip-Up Cone Lock X-Large-----	105FP17	\$170
G <input type="checkbox"/> One Piece Footboard (24"-32")-----	105FP20	\$195

17 FOOTREST /LEGREST ACCESSORIES

<input type="checkbox"/> Heel Loops -----	105FA01	Standard
<input type="checkbox"/> Omit Heel Loops -----	105FA02	less net \$7
<input type="checkbox"/> Padded Calf Strap -----	105FA03	\$70
H <input type="checkbox"/> Swingaway Calf Pad (pair) -----	105FA04	\$185
<input type="checkbox"/> Foam Guard on Elevating Legrest (pair) -----	105FA05	\$70
<input type="checkbox"/> Foam Pad One Piece Footboard -----	105FA06	\$80
<input type="checkbox"/> Foam Pad on Footplate (pair) -----	105FA07	\$70

FOOTNOTES:

- F. Includes Swingaway Calf Pads
- G. Not available with 7" or 8" casters when footboard is facing back on 90 degree Swingaway
- H. Swingaway Calf Pad can only be assembled onto Swingaway front rigging.
- I. Only available with Stroller Bar push option.
- J. Not available with Flip Back Arms

19 WHEEL LOCKS

<input type="checkbox"/> Push-to-Lock Wheel Lock -----	105WL1	Standard
<input type="checkbox"/> Pull-to-Lock Wheel Lock -----	105WL2	\$40
<input type="checkbox"/> Wheel Lock Extension -----	105WL3	No Charge

20 HEADREST HARDWARE OPTIONS

<input type="checkbox"/> Fixed Headrest Hardware -----	105HE01	Standard
<input type="checkbox"/> Fixed Headrest Hardware Ass'y (Plus 4")-----	105HE02	\$20
<input type="checkbox"/> Single Multi-Axial Headrest Hardware -----	105HE03	\$158
<input type="checkbox"/> Single Multi-Axial Headrest Hardware (Plus 4") -----	105HE04	\$183
<input type="checkbox"/> Dual Multi-Axial Headrest hardware -----	105HE05	\$195
<input type="checkbox"/> Dual Multi-Axial Headrest hardware (Plus 4") -----	105HE06	\$220
<input type="checkbox"/> Links Multi Adjustable Headrest -----	105HE10	\$200
<input type="checkbox"/> Links Multi Adjustable Curved Headrest -----	105HE11	\$240
<input type="checkbox"/> Additional Link for Headrest -----	105HE12	\$50
<input type="checkbox"/> 2 Additional Links for Headrest -----	105HE13	\$100
<input type="checkbox"/> 3 Additional Links for Headrest -----	105HE14	\$150
<input type="checkbox"/> UHRBA -----	105HEBA	\$61

② HEADREST FOAM PAD OPTIONS

<input type="checkbox"/> Small Headrest Foam Pad -----	105HF1	No Charge
<input type="checkbox"/> Small Headrest with Wings Foam Pad -----	105HF2	\$90
<input type="checkbox"/> Medium Headrest Foam Pad -----	105HF7	No Charge
<input type="checkbox"/> Medium Headrest with Wings Foam Pad -----	105HF8	\$90
<input type="checkbox"/> Large Headrest Foam Pad -----	105HF3	Standard
<input type="checkbox"/> Large Headrest with Wings Foam Pad -----	105HF4	\$95
<input type="checkbox"/> Omit Headrest -----	105HE09	less net \$10

② ANTI TIPPERS

<input type="checkbox"/> Fixed Anti-Tipper -----	105MA22	No Charge
--	---------	-----------

③ MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Auto Seat Positioning Strap -----	105MA03	\$52
<input type="checkbox"/> Airline Seat Positioning Strap -----	105MA04	\$76
<input type="checkbox"/> Velcro Release Positioning Strap -----	105MA05	\$48
<input type="checkbox"/> Foam Pad on Arm Guard (pair) -----	105MA06	\$140
<input type="checkbox"/> Padded Lap Tray 23"-24" -----	105MA51	\$385
<input type="checkbox"/> Padded Lap Tray 25"-26" -----	105MA52	\$385
<input type="checkbox"/> Padded Lap Tray 27"-28" -----	105MA53	\$385
<input type="checkbox"/> Padded Lap Tray 29"-30" -----	105MA54	\$385
<input type="checkbox"/> Padded Lap Tray 31"-32" -----	105MA55	\$385
<input type="checkbox"/> Half Padded Lap Tray, left side-----	105MA14	\$160
<input type="checkbox"/> Half Padded Lap Tray, right side-----	105MA15	\$160
<input type="checkbox"/> Elbow Stop, left -----	105MA16	\$94
<input type="checkbox"/> Elbow Stop, right -----	105MA17	\$94
ⁱ <input type="checkbox"/> Shoulder Lateral Pad (each) -----	105MA18	\$90
^j <input type="checkbox"/> Adjustable Side Pad -----	105MA19	\$157
<input type="checkbox"/> IV Pole, telescopic-----	105MA20	\$168
<input type="checkbox"/> Oxygen Tank Holder-----	105MA21	\$173
<input type="checkbox"/> Quick Release Axle -----	105AL1	\$110

WARRANTY

Future Mobility Healthcare Inc. warrants the Orion II 500 to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.