



<b>Account Information</b>	Order Date: _____
Account #: _____	P.O.# _____
Business Name: _____	Phone: _____
Ship to Address: _____	City: _____
Fax: _____	Province: _____ Postal Code: _____
Contact: _____	Tagged For: _____

<b>① FRAME TYPE</b>		
<input checked="" type="checkbox"/> Orion II - WMA5FM005 (User Weight Capacity: 350lbs) -----	101	<b>\$2,700</b>
<b>② FRAME FINISH</b>		
<input checked="" type="checkbox"/> Black -----	101C01	No Charge
<b>③ SIZE CONFIGURATION</b>		
<input type="checkbox"/> 16w x 16D, 22"rear wheels and 6" casters (STF 16")	101QS1	No Charge
<input type="checkbox"/> 16w x 18D, 22"rear wheels and 6" casters (STF 16")	101QS2	No Charge
<input type="checkbox"/> 18w x 16D, 22"rear wheels and 6" casters (STF 16")	101QS3	No Charge
<input type="checkbox"/> 18w x 18D, 20"rear wheels and 5" casters (STF 15")	101QS4	No Charge
<input type="checkbox"/> 18w x 18D, 22"rear wheels and 6" casters (STF 16")	101QS5	No Charge
<input type="checkbox"/> 18w x 20D, 22"rear wheels and 6" casters (STF 16")	101QS6	No Charge
<input type="checkbox"/> 20w x 18D, 22"rear wheels and 6" casters (STF 16")	101QS7	No Charge
<b>④ SEAT DEPTH</b>		
<input checked="" type="checkbox"/> EZ Tilt Option -----	101AL4	No Charge
<b>⑤ REAR WHEEL RIM</b>		
<input checked="" type="checkbox"/> Mag Wheels -----	101RW09	No Charge
<b>⑥ REAR WHEEL TIRE</b>		
<input checked="" type="checkbox"/> Airless (Solid Rubber) -----	101RT1	Standard
<b>⑦ HANDRIM OPTIONS</b>		
<input checked="" type="checkbox"/> Chrome -----	101HR1	Standard
<b>⑧ ARM TYPE</b>		
<input checked="" type="checkbox"/> T-Style Adjustable Height -----	101AR04	Standard
<b>⑨ ARM PAD</b>		
<input checked="" type="checkbox"/> Contour Gel 17" (pair)-----	101AP3	Standard
<b>⑩ PUSH OPTIONS</b>		
<input checked="" type="checkbox"/> Canes & Stroller Bars Fixed High -----	101SP2	No Charge
<b>⑪ TILT TRIGGER LOCATION</b>		
<input checked="" type="checkbox"/> Trigger Location - Stroller Bar-----	101TH2	No Charge
<b>⑫ BACK OPTION</b>		
<input checked="" type="checkbox"/> Prism Basic Back -----	101BK01	\$150
<b>⑬ SEAT OPTION</b>		
<input checked="" type="checkbox"/> Prism Supreme Cushion -----	101SC07	\$375
<b>⑭ FRONT RIGGING</b>		
<input checked="" type="checkbox"/> 70 Degree Swingaway Upgraded -----	101FR04	\$215
<b>⑮ FOOTPLATE TYPE</b>		
<input checked="" type="checkbox"/> Aluminum Angle Adj. Flip Up Footplates Medium-----	101FP11	\$120
<b>⑯ FOOTREST /LEGREST ACCESSORIES</b>		
<input checked="" type="checkbox"/> Heel Loops -----	101FA01	Standard
<b>⑰ WHEEL LOCKS</b>		
<input checked="" type="checkbox"/> Push-to-Lock Wheel Lock -----	101WL1	Standard
<input type="checkbox"/> Wheel Lock Extension -----	101WL3	No Charge
<b>⑱ HEADREST HARDWARE OPTIONS</b>		
<input checked="" type="checkbox"/> Fixed Headrest Hardware -----	101HE01	No Charge
<input type="checkbox"/> Dual Multi Axial Headrest Hardware -----	101HE05	\$195
<b>⑲ HEADREST FOAM PAD OPTIONS</b>		
<input checked="" type="checkbox"/> Large Headrest Foam Pad -----	101HF3	No Charge
<input type="checkbox"/> Pro-Fit Medium Headrest -----	101HF13	\$105
<b>⑳ ANTI TIP</b>		
<input checked="" type="checkbox"/> Quick Flip Anti-Tippers -----	101MA23	Standard

<b>Standard Features:</b>
Tilt Function (45 Degrees)
Seat Width 18"
Seat Depth 18"
Stroller Bar Fixed Dual Bend
T-Style Height Adjustable Full Length Gel
Seat to Floor 16"
Rear Wheels 22"
Casters 6"
Multi Adjustable Axel Plate
70 Degree Front Rigging w/aluminum flip up footplate
Heel Loops
Push to Lock Wheel Locks
Aluminum Seat Pan
Prism Ideal Cushion
Adjustable Orion Back
Two way Padded, Contoured Headrest
Anti-tippers
Colour - Black



ORION II FRONT STF HEIGHT AVAILABILITY				
REAR WHEEL	4" CASTER			
16"	13.25"			
20"	13.25"			
REAR WHEEL	5" CASTER			
16"	14.5"			
20"	14.5"	15"	16"	
22"	15.5"			
REAR WHEEL	6" CASTER			
12"	15.5"			
20"	16"			
22"	16"			
24"	16"	16.5"		
REAR WHEEL	7" CASTER			
22"	16.5"			
24"	16.5"	17"		
REAR WHEEL	8" CASTER			
16"	18"			
24"	18"	19"	20"	