

# REHAB REVOLUTION

## REHAB COMMODE SERIES

 PROUDLY MADE IN **CANADA**

### FEATURES

- Foam Skinned or Mesh Back
- Quick Flip Arm Rests
- Quick Release Back Clips
- Easy Removable Commode Pan



### FEATURES

- Foam Skinned or Mesh Back
- Removable Back Frame
- Self-Propelled Option
- Easy Removable Commode Pan



### FEATURES

- Foam Skinned or Mesh Back
- Slide-Flip Arm Rests
- Quick Release Back Clips
- Easy Removable Commode Pan



### FEATURES

- Foam Skinned or Mesh Back
- Quick Flip Rigid Armrests
- Quick Release Back Clips
- Easy Removable Commode Pan
- 3-38° Of Tilt



### FEATURES

- Foam Skinned or Mesh Back
- Quick Flip Arm Rests
- Self-Propelled Option
- Full Side Access for Wiping
- Easy Removable Commode Pan



# AQUA

## COMMUNE - TL

### FEATURES

- Foam -Filled Rigid Cushion
- Slide-Flip Armrests
- Comfort Cushion
- Easy Remove Potty
- 3-38° Of Tilt



## NEW FOLDING COMMODE

### FC

Product Weight	70 (LBS)
Weight Capacity	350 (LBS)
Seat Width	14,16,18,20,22 (Inch)
Seat Depth	18"
Overall Width	Chair Width +7"
Overall Length	36-38"
Seat To Floor Height	17-20"
Back Height	19-22"
Caster Options	5"
Wheel Size	20,22,24 (Inch)
Wheel Type	Urethane
Arm Type	Fixed Cantilever Flip Up
Arm Height	8"
Width Between Armrests	18,20,22 (Inch)
Front Rigging	70° Removable 70° Swing-Away 50° Removable 50° Swing-Away
Degree of Tilt	3-33°
Frame Material	Stainless Steel
Frame Finish	Powder Coat

## FULL SPECIFICATION COMMODE LIST

	TL	TLF	NT*	NTF*	NT600
Product Weight	70 (LBS)	57 (LBS)	38 (LBS)	38 (LBS)	50 (LBS)
Weight Capacity	350 (LBS)	350 (LBS)	350 (LBS)	350 (LBS)	600 (LBS)
Seat Width	14",16",18",20",22"	14",18",20",22"	16",18",20",22"	18",20",22"	24",26"
Seat Depth	18"	18"	18"	18"	18"
Overall Width	Chair Width +6"	Chair Width +6"	Chair Width +6"	Chair Width +6"	Chair Width +6"
Overall Length	39"	38"	36-38"	36-38"	39"
Seat To Floor Height	20-23"	20-23"	17-20" (Self Prop) 20-23" (Auto Prop)	17-20" (Self Prop) 20-23" (Auto Prop)	20-23"
Back Height	17-20"	17-20"	19-22"	19-22"	19-22"
Caster Options	5"	5"	5"	5"	5"
Wheel Size (Self Propelled)	N/A	N/A	20,22,24 (Inch)	20,22,24 (Inch)	N/A
Wheel Type	Urethane	Urethane	Urethane*	Urethane*	Urethane
Arm Type	Slide Flip-Up Arms	Fixed Cantilever	Fixed Cantilever Fixed Flip Up	Fixed Cantilever Fixed Flip Up	Fixed Cantilever Fixed Flip Up
Arm Height	7"	8"	8"	8"	9"
Width Between	18,20,22"	16,18,20,22"	18,20,22"	18,20,22"	21.5",23.5"
Front Rigging	50° Removable 60° Swing Away	50° Removable 60° Swing Away	70° & 60° Swing-Away	70° & 60° Swing-Away	70° & 60° Swing-Away
Degree of Tilt	3-38°	3-38°	N/A	N/A	N/A
Frame Material	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel	Stainless Steel
Frame Finish	Powder Coat	Powder Coat	Powder Coat	Powder Coat	Powder Coat

### \* self-propelled option notes:

- Self Propelled Models add +7" To Overall Width Instead of +6"
- Pneumatic or Urethane Wheels Available with Self-Propelled Models