



CAPELLA DYNAMIC TILT

OPTIONS: 20 degrees of tilt (optional 45 degrees of tilt with Capella 45), 21 degrees of recline (optional), growable frame – Seat width can increase from 16" to 20" or 20" to 24". The seat's depth can increase by as much as 1" without the need to change frame members.

COLOR OPTIONS



LOWEST SEAT-TO-FLOOR-HEIGHT ON THE MARKET




40 LBS. WEIGHT **20° & 45° TILT DEGREE** **350 LBS. CAPACITY**

This Dynamic Tilt is designed with the lowest Seat-To-Floor Height on the market 11.75". It's Ideal for foot propellers



Scan me





11.75" SEAT-TO-FLOOR HEIGHT

QUICK RELEASE

Allows for quicker and easier changing or removing of the wheelchair tires. This option is of value to those who use the folding chair for transport and want to reduce weight of the chair when loading it into a vehicle.

ADJUSTABLE CASTER HEIGHT

There are different type of caster height-to floor options available. Select the right height according to the requirement.



ABILITY TO GROW

Frame grows as required without replacing any parts up to 4" of width growth. Seat depth is grow able up to 2" with each seat rail at no extra cost.

GAS CYLINDER

Single cylinder linkage mounts onto dual slide rails for ultra rigid frame and simplified design. It can be found located underneath the chair, allows for a gentle, slow tilt for the client lockable at infinite angles




SPECIFICATIONS:

CAPELLA DYNAMIC TILT



SEAT TO FLOOR

11 3/4" - 19"
(13" - 19" ON CAPELLA 45°)




REAR WHEEL

12" - 16", 20", 22" OR 24"
COMPOSITE URETHANE OR
SPOKED WHEELS



SEAT WIDTH

16" - 24"



SEAT DEPTH

15" - 20"



OVERALL WIDTH

11 3/4" - 19"
(13" - 19" ON CAPELLA 45)



CASTER SIZE

4" (ONLY ON CAPELLA)
5", 6", 7", 8" URETHANE



FRONT RIGGING

SWINGAWAY
DUAL SWINGAWAY
OR DUAL ELEVATING LEGRESTS



PRODUCT WEIGHT

66 LBS (30KG)



ADJUSTABLE BACK ANGLE

-7°, 0°, 7°, 14°, 21°




BACK STYLE

ORON BACK



ARM STYLE

ADJ. HEIGHT T-STYLE
ADJ. HEIGHT PIN RELEASE T-STYLE
ADJ. HEIGHT FLIP BACK



FRAME MATERIAL

ALUMINUM



WEIGHT CAPACITY

CAPELLA 20° - 250 LBS (113 KG)
CAPELLA 45° - 250 LBS (113 KG)



ADP
Ontario

CODE

WM0000559

WARRANTY:

The Capella, Capella45, Orion II, Orion III and the Super Chair for Kids is warranted to be free from defects in materials and workmanship for a period of one (1) year and the Orion II 500 for period of one (1) year from the date of purchase. The frame has a lifetime warranty and the upholstery for 90 days, upon normal usage by the original purchaser. Power tilt batteries are warranted for six (6) months while the control box, battery charger, actuator and hand control are warranted for one (1) year. If within this warranty period the product shall be proven to be defective, such product shall be repaired or replaced, at Future Mobility Healthcare Inc. discretion.



CANADA

6750 Professional Court, Mississauga,
Ontario, L4V 1X6



1-888-737-4011



www.futuremobility.ca



USA



1 Buffalo River Place, Buffalo,
NY, 14210



1-855-938-8873



www.futuremobility.com

