


Prism Ideal

The Prism Ideal Cushion ™ was designed for users that require mild to moderate postural support and offers gentle contours, helping to promote correct seating posture. A seat cushion with 3 distinct layers: a top layer of soft, comfortable poly foam; a central layer of pressure relieving visco-elastic foam; and a base layer of medium density polyurethane support foam, with a contour cut to provide positioning. We have created an economical and lightweight cushion with comfort, positioning and pressure management.



Stock Sizes - Same Day Shipment

ADP CODE: SESFM0050					
Part Number	Size	Part Number	r Size		
IC1616	16" x 16"	IC1818	18" x 18"		
IC1618	16" x 18"	IC2016	20" x 16"		
IC1816	18" x 16"	IC2018	20" x 18"		

Non-Stock Sizes - Shipment within 3 business days

ADP CODE: SESFM0050						
Part Number	Size	Part Number	Size			
IC1620	16" x 20"	IC1720	17" x 20"			
IC1716	17" x 16"	IC1820	18" x 20"			
IC1718	17" x 18"	IC2020	20" x 20"			

Product Highlights

- · Comfort Cushion for light everyday use
- Standard water resistant outer cover, four way stretch
- Waterfall front edge for foot propelling comfort
- Core layer of visco-elastic foam
- · Lightweight cushion
- Low profile, high-density foam base enables

- proper positioning
- · Three distinct layers of foam


Product Weight: 2 lbs User Weight: 250 lbs


Warranty: Two (2) years limited on cushion and 90 days on the cover.

Best Used For:


Adv. Weight Distribution
Pediatric Capability
Self-Adjustment
Mild Positioning Needs
Low Skin Breakdown Risk
Short Stint Use

Pressure Reading Scale in mmHG





Pressure Map	Prism Ideal
Weight	125 Lbs/57 Kg
Duration	1 hr
Average Pressure	32 mmHg



Pressure Map	Prism Ideal
Weight	200 Lbs/91 Kg
Duration	1 hr
Average Pressure	33 mmHg

CUSHION THICKNESS

A- Height at adductor = 3"

B- Height at abductor = 3"

C- Height at rear IT = 21/4"

D- Leg Trough = **2** 1/4"



CANADA

6750 Professional Court, Mississauga, L4V1X6

1-888-737-4011

www.futuremobility.ca

USA

1 Buffalo River Place, Buffalo, NY 14210

1-855-938-8873

www.futuremobility.com