




Supreme HD

The Prism Supreme Heavy Duty™ Cushion was designed for users weighing up to 500 lbs and requiring high positioning support. The anatomically designed, firm polyurethane layer, with its aggressive leg troughs and abductor contouring provides superior pelvic stability and leg positioning. A base layer of pressure relieving visco-elastic foam under the IT's (Ischial Tuberosity), allows for excellent immersion for bony prominences.



Custom Sizes - Shipment within 7 business days

ADP CODE: SESFM0125

Part Number	Size
SCHDCUS	Custom

Stock Sizes – Same Day Shipment

ADP CODE: SESFM0125

Part Number	Size	Part Number	Size
SCHD2416	24" x 16"	SCHD2820	28" x 20"
SCHD2418	24" x 18"	SCHD3016	30" x 16"
SCHD2420	24" x 20"	SCHD3018	30" x 18"
SCHD2616	26" x 16"	SCHD3020	30" x 20"
SCHD2618	26" x 18"	SCHD3216	32" x 16"
SCHD2620	26" x 20"	SCHD3218	32" x 18"
SCHD2816	28" x 16"	SCHD3220	32" x 20"
SCHD2818	28" x 18"		

Product Highlights

- High density, contour cut layered foam for a consistent, long lasting product
- Cane cut in rear to allow for correct placement in wheelchair and proper seat depth
- Radius cut to match shape of back upholstery
- Layer of visco foam along the bottom for further pressure relief under IT's (Ischial Tuberosity)
- Supreme pressure relief

- 4-way stretch, water resistant inner liner
- 4-way stretch, breathable outer cover

Product Weight: 7 lbs

User Weight: 500 lbs


Warranty: Two (2) years limited on cushion and 90 days on the cover.

Best Used For:
 Adv. Weight Distribution
 500 LBS +
 Extended Length Use
 Adv. Positioning Needs
 Adv. Immersion & Resilience
 Adv. Pelvic Stability

Non-Stock Sizes
 Shipment within 3 business days

ADP CODE: SESFM0125

Part Number	Size	Part Number	Size
SCHD2422	24" x 22"	SCHD2824	28" x 24"
SCHD2424	24" x 24"	SCHD2916	29" x 16"
SCHD2516	25" x 16"	SCHD2918	29" x 18"
SCHD2518	25" x 18"	SCHD2920	29" x 20"
SCHD2520	25" x 20"	SCHD2922	29" x 22"
SCHD2522	25" x 22"	SCHD2924	29" x 24"
SCHD2524	25" x 24"	SCHD3022	30" x 22"
SCHD2622	26" x 22"	SCHD3024	30" x 24"
SCHD2624	26" x 24"	SCHD3116	31" x 16"
SCHD2716	27" x 16"	SCHD3118	31" x 18"
SCHD2718	27" x 18"	SCHD3120	31" x 20"
SCHD2720	27" x 20"	SCHD3122	31" x 22"
SCHD2722	27" x 22"	SCHD3124	31" x 24"
SCHD2724	27" x 24"	SCHD3222	32" x 22"
SCHD2822	28" x 22"	SCHD3224	32" x 24"



CUSHION THICKNESS


A- Height at adductor = 6 1/2"

B- Height at abductor = 5"

C- Height at rear IT = 4 1/4"

D- Leg Trough = 4"

E- Height at rear trochanter = 5 1/4"



FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga, L4V1X6

1-888-737-4011

www.futuremobility.ca

USA

1 Buffalo River Place, Buffalo, NY 14210

1-855-938-8873

www.futuremobility.com