



# FUTURE MOBILITY HEALTHCARE INC.

DIVISION OF ABF MFG GROUP INC.  
WWW.FUTUREMOBILITY.CA



TEL: 905-671-1661  
FAX: 905-671-3377  
orderdesk@futuremobility.ca

### Account Information

Account #: \_\_\_\_\_  
Business Name: \_\_\_\_\_  
Ship to Address: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Tag For: \_\_\_\_\_

Order Date: \_\_\_\_\_  
P.O.# \_\_\_\_\_  
Phone: \_\_\_\_\_  
City: \_\_\_\_\_  
Province: \_\_\_\_\_ Postal Code: \_\_\_\_\_

### 1 FRAME TYPE

- Attendant Propelled Stainless Steel (white) 107C19 \$1,345
- Self Propelled Stainless Steel (white) 107C19S \$1,895

### 2 SEAT WIDTH

- 17" 107SW17 No Charge
- 20" 107SW20 \$150
- 22" 107SW22 \$200

### 3 SEAT DEPTH

- 17" 107SD17 No Charge

### 4 SEAT-TO-FLOOR HEIGHT (ON ATTENDANT PROPELLED COMMODE)

- 16.5" 107SF09 No Charge
- 17.5" 107SF11 Standard
- 18.5" 107SF13 No Charge
- 19.5" 107SF15 No Charge

\* Seat to Floor measure without seatpad. (Thickness 2.75")

### 5 CASTER SIZE

- 5" x 2" Polypropylene 107CW03 Standard

### 6 CASTER TYPE

- Total Lock 107CT01 Standard

### 7 WHEEL SIZE (FOR SELF PROPELLED COMMODE)

- 20" Composite Urethane Tires 107RW03 No Charge  
Seat to Floor with 20" wheels is 16" or 18"
- 22" Composite Urethane Tires 107RW04 Standard  
Seat to Floor with 22" wheels is 17" or 19"
- 24" Composite Urethane Tires 107RW05 No Charge  
Seat to Floor with 24" wheels is 18" or 20"

### 8 HANDRIM OPTIONS (FOR SELF PROPELLED COMMODE)

- Chrome 107HR1 Standard
- Plastic Coated 107HR2 \$160
- Omit Handrims 107HR3 Less Net \$15


### 9 ARM TYPE

- Cantilever Arm 107AR29 Standard
- Flip Back Fixed Height 107AR23 \$250

### 10 PUSH OPTIONS

- Stroller Handles 107SH1 Standard

Commode Type	STF Front	Caster Hole	Rear Wheel Hole	Rear Wheel	STF Rear
Attendant Prop.	16.5"	D	N/A	N/A	16.5"
Attendant Prop.	17.5"	C	N/A	N/A	17.5"
Attendant Prop.	18.5"	B	N/A	N/A	18.5"
Attendant Prop.	19.5"	A	N/A	N/A	19.5"
Self Propelled	16.25"	D	A	20"	15.75"
Self Propelled	17.25"	C	A	22"	16.5"
Self Propelled	18.25"	B	B	20"	17.75"
Self Propelled	18.25"	B	A	24"	17.5"
Self Propelled	19.25"	A	B	22"	18.5"
Self Propelled	19.5"	A	B	24"	19.25"



### Standard Features:

- Seat Width 17"
- Seat Depth 17"
- Adjustable Height Back
- Slip on Back Support
- Seat to Floor Height 16.5" to 19.5" adj. without tools on Attendant Propelled w/4 casters
- Casters 5" - Total Lock - Stainless Steel
- Cantilever Fixed Height Arm Supports
- Removable Seat ( Thickness 2.75" )
- Removable Pan Support
- Removable Pan with lid
- Finish - White Powder Coat

### 11 BACK OPTIONS

- Slip on Mesh Back Support 107BK46 Standard
- Foam Skinned Back Support 107BK47 No Charge

### 12 SEAT OPTIONS

- Removable Seat 107SO1 Standard
- Removable Pan Support 107SO2 Standard
- Removable Pan with Lid 107SO3 Standard

### 13 FRONT RIGGING

- 50° Swingaway Stainless Steel (white) 107FR13 \$195
- 70° Swingaway Stainless Steel (white) 107FR22 \$210

### 14 FOOTPLATE TYPE

- Composite Flip-Up 107FP02 No Charge
- Composite Adj. Angle Flip-Up Spring Return 107FP04 \$150

### 15 FOOTREST/LEGREST ACCESSORIES

- Heel Loops 107FA01 No Charge
- Omit Heel Loops 107FA02 Less Net \$7

### 16 WHEEL LOCKS (FOR SELF PROPELLED COMMODE)

- Push-to-Lock 107WL1 Standard
- Wheel Lock Extension 107WL3 \$35

### 17 ANTI TIP

- Rear Anti-Tippers 107MA22 Standard

### 18 MISCELLANEOUS ACCESSORIES (FOR SELF PROPELLED COMMODE)

- Quick Release Axles 107AL1 \$110

### WARRANTY:

Future Mobility Healthcare warrants its Aqua Commode NT to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the two (2) years of the original purchaser/user and the upholstery for 90 days, upon normal usage by original purchaser. If within this warranty period the product shall be proven to be defective, such product shall be repaired or replaced, at the discretion of Future Mobility Healthcare Inc

### NOTES: