



FUTURE MOBILITY HEALTHCARE INC.

DIVISION OF ABF MFG GROUP INC.
WWW.FUTUREMOBILITY.CA



NORTHERN LITE

by FUTURE MOBILITY

TEL: 905-671-1661
FAX: 905-671-3377
orderdesk@futuremobility.ca

Account Information

Account #: _____
Business Name: _____
Ship to Address: _____
Contact: _____
Tag For: _____

Order Date: _____
P.O.# _____
Phone: _____
City: _____
Province: _____ Postal Code: _____

1 FRAME TYPE

Northern Lite - WMA3FM005 108 \$2,350
^A Front seat height should be larger than rear seat height

2 FRAME FINISH

Silver Lining 108C02 Standard
 Wet Black 108C01 No Charge

3 SEAT WIDTH

16" 108SW10 No Charge
 18" 108SW18 No Charge
 20" 108SW20 No Charge

4 SEAT DEPTH

16" 108SD16 No Charge
 18" 108SD18 No Charge
 20" 108SD20 No Charge

5 SEAT-TO-FLOOR HEIGHT AND CASTER SIZE

13.75" (4" Caster) 108STFN1 No Charge
 14.5" (4" Caster) 108STFN2 No Charge
 15.75" (4" Caster) 108STFN3 No Charge
 16.5" (4" Caster) 108STFN4 No Charge
 16" (6" Caster) 108STFN5 No Charge
 16.75" (6" Caster) 108STFN6 Standard
 17.5" (6" Caster) 108STFN7 No Charge
 17.75" (8" Caster) 108STFN8 No Charge
 18.5" (8" Caster) 108STFN9 No Charge
 19" (8" Caster) 108STFN10 No Charge

6 REAR SEAT-TO-FLOOR HEIGHT AND WHEEL SIZE

^B 13" (20" Wheels) 108STFN11 No Charge
 13.75" (20" Wheels) 108STFN12 No Charge
^B 14" (22" Wheels) 108STFN13 No Charge
 14.5" (20" Wheels) 108STFN14 No Charge
 14.75" (22" Wheels) 108STFN15 No Charge
 15.25" (20" Wheels) 108STFN16 No Charge
 15.5" (22" Wheels) 108STFN17 No Charge
 15.75" (24" Wheels) 108STFN18 No Charge
 16.25" (22" Wheels) 108STFN19 No Charge
 16.5" (24" Wheels) 108STFN20 Standard
 17.25" (24" Wheels) 108STFN21 No Charge
 18" (24" Wheels) 108STFN22 No Charge
 18.75" (24" Wheels) 108STFN23 No Charge

^B Can only have 0° camber angle on rear wheels at these heights

7 CASTER OPTIONS

4" x 1.25" Urethane 108CW01 No Charge
 5" x 1.25" Urethane 108CW05 No Charge
 6" x 1.25" Urethane 108CW04 No Charge
 8" x 1.25" Urethane 108CW07 No Charge

8 REAR WHEEL TIRE

Airless (Solid Rubber) 108RT1 Standard
 Pneumatic 108RT2 No Charge

Standard Features:

- Lightweight Folding Frame w/Swingaway Front End
- Elliptical Heavyduty Double Crossbrace Tubing
- Seat Width (16"-20")
- Seat Depth (16"-18")
- Front Seat-To-Floor Height (14"-19")
- Rear Wheel (20", 22" or 24") Airless
- Caster Options (4", 6" or 8")
- Camber Adj. Rear Wheels (0°, 3°, 5°)
- Seat Upholstery
- Adjustable Tension Back Upholstery
- Push Canes
- Composite Flip-Up Footplates
- Heel Loops
- Wheel Locks (Push or Pull to Lock)
- Full Length Adj. Height T-Arm (Vinyl Pad)
- Anti tippers

9 REAR WHEEL OPTIONS

<input type="checkbox"/> 20" Spoke Wheels	108RW06	No Charge
<input type="checkbox"/> 22" Spoke Wheels	108RW07	No Charge
<input type="checkbox"/> 24" Spoke Wheels	108RW08	No Charge
<input type="checkbox"/> 20" Mag Wheels	108RW09	No Charge
<input type="checkbox"/> 22" Mag Wheels	108RW10	No Charge
<input type="checkbox"/> 24" Mag Wheels	108RW11	Standard

10 HANDRIM OPTIONS

<input type="checkbox"/> Chrome	108HR1	Standard
<input type="checkbox"/> Plastic Coated	108HR2	\$110
<input type="checkbox"/> Omit Handrims	108HR3	Less Net \$15


11 ARM TYPE

<input type="checkbox"/> T-Style Adjustable Height Pin Release - Adult	108AR06	No Charge
<input type="checkbox"/> T-Style Adjustable Height Pin Release - Youth	108AR05	No Charge
^C <input type="checkbox"/> Flip Back Adjustable Height	108AR09	No Charge

^C Not available with Fold Down Cane

Armrest Matrix for -T Style Pin Release Arm

Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range
Adult	6.75	7.50	9.50" - 14.50"
Youth	4.50	5.0	8.75" - 13.75"



Arm Pad Length L

Desk	Full	Curved
10.25"	14.0"	17.0"

12 ARM PAD

^D <input type="checkbox"/> Desk length Vinyl 10" (pair)	108AP1	No Charge
<input type="checkbox"/> Full length Vinyl 14" (pair)	108AP2	No Charge
^E <input type="checkbox"/> Contour Gel 17" (pair)	108AP3	Standard

^D Not available with Flip Back Arm

^E Only Gel Arm Pad on Flip Back Arm when Narrow Chair Width Option used

13 BACK TYPE

^F <input type="checkbox"/> Back 10° Bend w/Push Handle Low	108PC3	Standard
^G <input type="checkbox"/> Back 10° Bend w/Push Handle High	108PC4	No Charge

^F Top of back cane is set at 16" high from the top of the seat rail

^G Top of back cane is set at 20" high from the top of the seat rail

14 BACK OPTIONS

<input type="checkbox"/> Slip-on Padded Back Upholstery	108BK48	Standard
<input type="checkbox"/> Adjustable Tension Back	108BK45	\$150
<input type="checkbox"/> Prism Basic Back	108BK01	\$150
<input type="checkbox"/> Prism Ultra Back	108BK06	\$250
<input type="checkbox"/> Prism Air Back	108BK11	\$400
<input type="checkbox"/> Prism TrueFitt Back	108BK16	\$650
<input type="checkbox"/> Prism Max Ultra Back	108BK21	\$300
<input type="checkbox"/> Prism Max Air Back	108BK26	\$450
<input type="checkbox"/> Prism Max TrueFitt Back	108BK31	\$700
<input type="checkbox"/> Lunar Back Mid Basic Lateral	108BK39	\$550
<input type="checkbox"/> Omit Back	108BK49	No Charge

W () x H ()

For other Backs, please fill in the required Backs Order Form and attach it to this order

15 BACK PINS (PRISM BACKS ONLY)

<input type="checkbox"/> Longer Pins for Back Hardware 1" Growth	108BP1	\$50
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth	108BP2	\$50

16 SEAT OPTIONS

<input type="checkbox"/> Black Nylon Seat Upholstery	108SC19	Standard
<input type="checkbox"/> Prism Ideal Cushion	108SC01	\$175
<input type="checkbox"/> Prism Elite Cushion	108SC04	\$250
<input type="checkbox"/> Prism Supreme Cushion	108SC07	\$375
<input type="checkbox"/> Prism Supreme II Plus Cushion	108SC10	\$375
<input type="checkbox"/> Prism Supreme Air Cushion	108SC13	\$500
<input type="checkbox"/> Omit Cushion	108SC20	No Charge

W () x D ()

Standard Size Cushion Depth: 16", 18" and 20"

For other Cushions, please fill in the required Cushions Order Form and attach it to this order

17 FRONT RIGGING

<input type="checkbox"/> 60° Swingaway Upgraded	108FR03	Standard
<input type="checkbox"/> 70° Swingaway Upgraded	108FR04	No Charge
<input type="checkbox"/> Swingaway Elevating Legrest (Pair)	108FR10	\$550

18 FOOTPLATE TYPE

<input type="checkbox"/> Composite Flip Up <small>Small</small>	108FP01	Standard
<input type="checkbox"/> Composite Flip Up <small>Large</small>	108FP02	No Charge
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Small</small>	108FP10	\$120
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Medium</small>	108FP11	\$120
<input type="checkbox"/> Aluminum Depth & Angle Adj. Flip-Up <small>Large</small>	108FP12	\$120

19 FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> Padded Calf Strap	108FA03	\$70
<input type="checkbox"/> Foam Pad on Elevating Legrest (Pair)	108FA05	\$70
<input type="checkbox"/> Heel Loops	108FA01	Standard
<input type="checkbox"/> Omit Heel Loops	108FA02	Less Net \$7

20 WHEEL LOCKS

^H <input type="checkbox"/> Push-to-Lock	108WL5	Standard
^I <input type="checkbox"/> Pull-to-Lock	108WL6	No Charge
<input type="checkbox"/> Wheel Lock Extension (Bolt on)	108WL7	No Charge

^H Standard Mount - Footrests swing further back to rear wheels^I Transfer Mount - Brakes levers are below seat level when locked**21 WHEEL LOCKS**

^J <input type="checkbox"/> Push-to-Lock	108CB1	Standard
<input type="checkbox"/> 3° Camber	108CB2	No Charge
<input type="checkbox"/> 5° Camber	108CB3	No Charge

^J Standard Wheelchairs are set at a 0° camber angle on rear wheels**20 ANTI TIP**

<input type="checkbox"/> Quick Flip Anti-Tippers	108MA23	No Charge
--	---------	-----------

21 MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Auto Seat Positioning Strap	108MA03	\$52
<input type="checkbox"/> Velcro Release Positioning Strap	108MA05	\$48

WARRANTY:

Future Mobility Healthcare warrants the Northern Lite wheelchair to be free from defects in materials and workmanship for a period of two (2) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.

NOTES: