



# FUTURE MOBILITY HEALTHCARE INC.

DIVISION OF ABF MFG GROUP INC.  
WWW.FUTUREMOBILITY.CA



TEL: 905-671-1661  
FAX: 905-671-3377  
orderdesk@futuremobility.ca

## Account Information

Account #: \_\_\_\_\_  
Business Name: \_\_\_\_\_  
Ship to Address: \_\_\_\_\_  
Contact: \_\_\_\_\_  
Tag For: \_\_\_\_\_

Order Date: \_\_\_\_\_  
P.O.# \_\_\_\_\_  
Phone: \_\_\_\_\_  
City: \_\_\_\_\_  
Province: \_\_\_\_\_ Postal Code: \_\_\_\_\_

## 1 FRAME TYPE

Galaxy Lite - WMA3FM010 103 \$2,350  
Seat angle is set at 3°

## 2 FRAME FINISH

Wet Black 103C01 No Charge  
 Silver Lining 103C02 No Charge  
**Custom Colours**  
 Cobalt Blue 103C03 \$150  
 Green Apple 103C04 \$150  
 Burnt Orange 103C05 \$150  
 Grape Purple 103C06 \$150  
 Fire Red 103C07 \$150  
 Yellow Ribbon 103C08 \$150  
 Custom Colours 103C17 \$150  
Allow 7-10 days for powder coating time for custom colors

## 3 SEAT WIDTH

\*  10" 103SW10 No Charge  
\*  11" 103SW11 No Charge  
\*  12" 103SW12 No Charge  
\*  13" 103SW13 No Charge  
\*  14" 103SW14 No Charge  
 15" 103SW15 No Charge  
 16" 103SW16 No Charge  
 17" 103SW17 No Charge  
 18" 103SW18 No Charge  
 19" 103SW19 No Charge  
 20" 103SW20 No Charge  
 22" 103SW22 \$250  
 24" 103SW24 \$250  
\*Must have Superchair Footplates on 10"-14" depth chair

## 4 SEAT DEPTH (EXTENDS FROM FRONT OF CHAIR)

\*  12" 103SD15 No Charge  
\*  13" 103SD16 No Charge  
\*  14" 103SD17 No Charge  
 15" 103SD18 No Charge  
 16" 103SD19 No Charge  
 17" 103SD20 \$110  
 18" 103SD21 \$250  
 19" 103SD22 \$250  
 20" 103SD23 \$250  
\*12"-14" not available with 22" or 24" rear wheels

Rear Wheel	4" Caster									
20"	13"	14"	14.5"	15"	16"					
22"				15.5"						
Rear Wheel	5" Caster									
20"				15"	15.5"	16.5"				
22"					16"	17"	17.5"			
Rear Wheel	6" Caster									
20"						17"	17.5"	18"		
22"							18"	18.5"		
24"							18"	18.5"		
Rear Wheel	7" Caster									
22"								19"	18.5"	

## Technical Specs:

- Lightweight Folding Frame
- Elliptical Heavy Duty Cross Brace Tubing
- Seat Width 10" - 20", 22" - 24"
- Seat Depth 12" - 20"
- 4", 5", 6", 7" Casters
- 20", 22", 24" Rear Wheels
- Angle Adjustable Push Canes
- Handrim Options (Chrome or Plastic Coated)
- Slip on Seat Upholstery
- Slip on Back Upholstery
- Arm rest height adjustable (1/2" increments)
- Push to Lock Wheel Locks
- Aluminum Flip-Up Footplate
- Heel Loops

## 5 SEAT-TO-FLOOR HEIGHT

Can not use a 8" caster with this wheelchair

STF	Wheels	Casters	Low Frame	STF Code	Price
<input type="checkbox"/> 13.0"	20"	4"	Low Frame	103STF45	\$225
Short fork is used to achieve these front STF heights					
<input type="checkbox"/> 14.0"	20"	4"	Low Frame	103STF46	No Charge
For front STF of 14" and above the long fork is used					
<input type="checkbox"/> 14.5"	20"	4"	Low Frame	103STF47	No Charge
<input type="checkbox"/> 15.0"	20"	4"	Low Frame	103STF48	No Charge
<input type="checkbox"/> 15.0"	20"	5"	Low Frame	103STF03	No Charge
<input type="checkbox"/> 15.5"	20"	5"	Low Frame	103STF04	No Charge
<input type="checkbox"/> 15.5"	22"	4"	Low Frame	103STF49	No Charge
<input type="checkbox"/> 16.0"	20"	4"	Low Frame	103STF50	No Charge
<input type="checkbox"/> 16.0"	22"	5"	Low Frame	103STF05	No Charge
<input type="checkbox"/> 16.5"	20"	5"	Low Frame	103STF51	No Charge
<input type="checkbox"/> 17.0"	20"	6"	STD Frame	103STF52	No Charge
<input type="checkbox"/> 17.0"	22"	5"	STD Frame	103STF53	No Charge
<input type="checkbox"/> 17.5"	20"	6"	STD Frame	103STF54	No Charge
<input type="checkbox"/> 17.5"	22"	5"	STD Frame	103STF55	No Charge
<input type="checkbox"/> 18.0"	20"	6"	STD Frame	103STF56	No Charge
<input type="checkbox"/> 18.0"	22"	6"	STD Frame	103STF57	No Charge
<input type="checkbox"/> 18.0"	24"	6"	STD Frame	103STF20	No Charge
<input type="checkbox"/> 18.5"	22"	6"	STD Frame	103STF58	No Charge
<input type="checkbox"/> 18.5"	24"	6"	STD Frame	103STF59	No Charge
<input type="checkbox"/> 19.0"	22"	7"	STD Frame	103STF60	No Charge
<input type="checkbox"/> 19.5"	22"	7"	STD Frame	103STF61	No Charge

Front STF	Rear Wheel	Caster Size	Fork Size	Fork Hole	Axle Plate Hole	Rear STF	Stem Size	Side Frame
13.00	20"	4	Short	A	Super Low	12.50"	Short	Low
14.00	20"	4	Long	A	1	13.50"	Short	Low
14.50	20"	4	Long	B	1	14.00"	Short	Low
15.00	20"	4	Long	C	2	14.25"	Short	Low
15.00	20"	5	Long	B	2	14.25"	Short	Low
15.50	20"	5	Long	C	2	14.75"	Short	Low
15.50	22"	4	Long	B	2	14.75"	Long	Low
16.00	20"	4	Long	C	3	15.25"	Long	Low
16.00	22"	5	Long	B	2	15.25"	Long	Low
16.50	20"	5	Long	C	3	16.00"	Long	Low
17.00	20"	6	Long	B	3	16.25"	Short	Standard
17.00	22"	5	Long	C	2	15.75"	Short	Standard
17.50	20"	6	Long	C	4	17.00"	Short	Standard
17.50	22"	5	Long	B	3	17.00"	Long	Standard
18.00	20"	6	Long	B	4	17.25"	Long	Standard
18.00	22"	6	Long	B	3	17.00"	Long	Standard
18.00	24"	6	Long	B	3	17.50"	Long	Standard
18.50	22"	6	Long	C	4	18.25"	Long	Standard
18.50	24"	6	Long	C	3	18.00"	Long	Standard
19.00	22"	7	Long	C	4	18.25"	Long	Standard
19.50	22"	7	Long	C	4	18.25"		Standard

6750 Professional Court, Mississauga, ON, L4V 1X6 Phone: +1 (888) 737-4011, (905) 671-1661 Fax: (866) 670-3377  
E-mail: orderdesk@futuremobility.ca Website: www.FutureMobility.ca

**6 REAR WHEEL RIM**

- Mag Wheels 103RW09 No Charge
- Spoke Wheels 103RW10 \$75

**7 REAR WHEEL TIRE**

- Airless (Solid Rubber) 103RT1 Standard
- Pneumatic 103RT2 No Charge

**8 CAMBER OPTIONS**

- 0° Camber 103CB1 Standard
- 3° Camber 103CB2 No Charge

**9 HANDRIM OPTIONS**

- Chrome 103HR1 Standard
- Plastic Coated 103HR2 \$110
- Omit Handrims 103HR3 less net \$15


**10 ARM TYPE**

- T-Style Adjustable Height Pin Release - Adult 103AR06 \$280
- T-Style Adjustable Height Pin Release - Youth 103AR05 \$300
- Swingaway Padded Tubular Arm 103AR01 \$250
- Flip Back Adjustable Height 103AR09 \$275

\*Not available with Fold Down Cane

**Armrest Matrix for -T Style Pin Release Arm**

Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range
Adult	6.75	7.50	9.50" - 14.50"
Youth	4.50	5.0	8.75" - 13.75"



**Arm Pad Length L**

Desk	Full	Curved
10.25"	14.0"	17.0"

**11 ARM PAD**

- Desk length Vinyl 10" (pair) 103AP1 No Charge  
Not available with Flip Back Arm
- Full length Vinyl 14" (pair) 103AP2 No Charge
- Contour Gel 17" (pair) 103AP3 Standard  
Only Gel Arm Pad on Flip Back Arm when Narrow Chair Width Option used

**12 BACK TYPE**

- Angle Adjustable 15° Bend Push Handle Low 103PC3 Standard
- Angle Adjustable 15° Bend Push Handle High 103PC4 \$90
- \*  Angle Adjustable 15° Bend Stroller Bar Low 103SP3 \$135
- \*  Angle Adjustable 15° Bend Stroller Bar High 103SP4 \$155

\*Stroller bar option only available for chair width 10"-15"

**13 BACK OPTIONS**

- Slip-on Padded Back Upholstery 103BK48 Standard
- Orion Back 103BK00 \$150
- Prism Basic Back 103BK01 \$150
- Prism Ultra Back 103BK06 \$250
- Prism Air Back 103BK11 \$400
- Prism TrueFitt Back 103BK16 \$650
- Prism Max Ultra Back 103BK21 \$300
- Prism Max Air Back 103BK26 \$450
- Prism Max TrueFitt Back 103BK31 \$700
- Lunar Back Mid Basic Lateral 103BK39 \$550
- Omit Back 103BK49 No Charge

W ( ) x H ( )

For other Backs, please fill in the required Backs Order Form and attach it to this order

**14 BACK PINS (PRISM BACKS ONLY)**

- Longer Pins for Back Hardware 1" Growth 103BP1 \$50
- Longer Pins for Back Hardware 2" Growth 103BP2 \$50

**NOTES:**

**WARRANTY:**

Future Mobility Healthcare warrants the Galaxy Lite wheelchair to be free from defects in materials and workmanship for a period of two (2) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.

**15 SEAT OPTIONS**

- Black Nylon Seat Upholstery 103SC19 Standard
- Prism Ideal Cushion 103SC01 \$100
- Prism Elite Cushion 103SC04 \$250
- Prism Supreme Cushion 103SC07 \$375
- Prism Supreme II Plus Cushion 103SC10 \$375
- Prism Supreme Air Cushion 103SC13 \$500
- Omit Cushion 103SC20 No Charge

W ( ) x D ( )

For other Cushions, please fill in the required Cushions Order Form and attach it to this order

**16 FRONT RIGGING**

- 70° Pin Style 103FR01 Standard
- 80° Pin Style 103FR02 \$135
- 60° Swingaway Upgraded 103FR03 \$215
- 70° Swingaway Upgraded 103FR04 \$215
- 90° Swingaway Upgraded 103FR05 \$230
- 60° Dual Swingaway Upgraded 103FR06 \$210
- 70° Dual Swingaway Upgraded 103FR07 \$210
- 80° Dual Swingaway Upgraded 103FR08 \$210
- 90° Dual Swingaway Upgraded 103FR09 \$210
- Dual Swingaway Elevating Legrest STD Frame Only 103FR25 \$550
- Swingaway Elevating Legrest (Pair) 113FR10 \$550
- Omit Front Rigging 103FR11 Less net \$20
- Omit Footrest Hanger 103FR12 Less net \$10

**17 FOOTPLATE TYPE**

- Composite Flip Up Small 103FP01 No Charge
- Composite Flip Up Large 103FP02 No Charge
- Aluminum Flip-Up Small 103FP05 Standard
- Aluminum Flip-Up Medium 103FP06 Standard
- Aluminum Flip-Up Large 103FP07 Standard
- Composite Angle Adj. Flip-Up Cone Lock Medium 18" 103FP09 \$120
- Composite Angle Adj. Flip-Up Cone Lock Large 103FP16 \$150
- Composite Angle Adj. Flip-Up Cone Lock X-Large 103FP17 \$170
- Aluminum Depth & Angle Adj. Flip-Up Small 103FP10 \$175
- Aluminum Depth & Angle Adj. Flip-Up Medium 103FP11 \$175
- Aluminum Depth & Angle Adj. Flip-Up Large 103FP12 \$175
- Composite Depth & Angle Adj. Small 103FP27 \$120
- Composite Depth & Angle Adj. Medium 103FP28 \$120
- Composite Depth & Angle Adj. Large 103FP29 \$120
- Composite Angle Adj. Flip-Up Spring Return 103FP04 \$150
- Superchair Kids Aluminum Flip-Up 103FP14 \$150  
Must have Superchair Footplates on 10"-14" depth chair
- One Piece Flip-Up Footboard 103FP19 \$160  
Not Available with 6" or 8" Casters when footboard is facing rear on 90° Swingaway

**18 FOOTREST/LEGREST ACCESSORIES**

- Padded Calf Strap 103FA03 \$70
- Foam Pad on Footplate (pair) 103FA07 \$70
- Foam Pad One Piece Footboard 103FA06 \$50
- Omit Heel Loops 103FA02 No Charge
- Heel Loops 103FA01 Standard

**19 WHEEL LOCKS**

- Push-to-Lock 103WL1 Standard
- Pull-to-Lock 103WL2 No Charge
- Wheel Lock Extension 103WL3 No Charge

**20 ANTI TIP**

- Quick Flip Anti-Tippers 103MA23 No Charge

**21 MISCELLANEOUS ACCESSORIES**

- Quick Release Axles 103AL1 No Charge
- Auto Seat Positioning Strap 103MA03 \$52
- Velcro Release Positioning Strap 103MA05 \$48
- Padded Lap Tray 15"- 16" 103MA09 \$285
- Padded Lap Tray 17"- 18" 103MA10 \$285
- Padded Lap Tray 19"- 20" 103MA11 \$285
- Half Padded Lap Tray, left side 103MA14 \$350
- Half Padded Lap Tray, right side 103MA15 \$350
- Arm Guard for Cantilever Arms 103MA50 \$110