



FUTURE

MOBILITY HEALTHCARE INC.

DIVISION OF ABF MFG GROUP INC.
WWW.FUTUREMOBILITY.CA



TEL: 905-671-1661
FAX: 905-671-3377
orderdesk@futuremobility.ca

Account Information

Account #: _____
Business Name: _____
Ship to Address: _____
Contact: _____
Tag For: _____

Order Date: _____
P.O.# _____
Phone: _____
City: _____
Province: _____ Postal Code: _____

1 FRAME TYPE

Superchair for kids - WMK5FM010 104 \$3,250

2 FRAME FINISH

Wet Black 103C01 No Charge
 Silver Lining 103C02 No Charge

Custom Colours

Cobalt Blue 103C03 No Charge
 Burnt Orange 103C05 No Charge
 Fire Red 103C07 No Charge
 Yellow Ribbon 103C08 No Charge
 Burgundy 104C11 No Charge
 Custom Colours 103C17 \$150

Allow 7-10 days for powder coating time for custom colors

3 SEAT WIDTH

10" 104SW10 No Charge
 11" 104SW11 No Charge
 12" 104SW12 No Charge
 13" 104SW13 No Charge
 14" 104SW14 No Charge
 15" 104SW15 No Charge
 16" 104SW16 No Charge

4 SEAT DEPTH

11" 104SD11 No Charge
 12" 104SD12 No Charge
 13" 104SD13 No Charge
 14" 104SD14 No Charge
 15" 104SD15 No Charge
 16" 104SD16 No Charge

5 SEAT-TO-FLOOR HEIGHT

17.5" (20" Wheels and 4" Casters) 104STF62 No Charge
 18" (20" Wheels and 4" Casters) 104STF63 No Charge
 18.5" (22" Wheels and 5" Casters) 104STF64 No Charge
 19" (22" Wheels and 7" Casters) 104STF60 No Charge
 19.5" (22" Wheels and 7" Casters) 104STF61 No Charge
 19.75" (22" Wheels and 8" Casters) 104STF65 No Charge
 18.5" (12" Wheels and 5" Casters) 104STF66 No Charge
 19" (12" Wheels and 7" Casters) 104STF67 No Charge
 19.5" (12" Wheels and 7" Casters) 104STF68 No Charge
 20" (12" Wheels and 8" Casters) 104STF69 No Charge

Front STF	Rear Wheel	Caster Size	Fork	Fork Hole	Axel Hole	Rear STF
18.5"	12"	4"	Long	C	A	17.25"
18.5"	12"	5"	Long	B	A	17.25"
19"	12"	7"	Short	C	A	17.5"
19.5"	12"	7"	Long	B	B	18.25"
20"	12"	8"	Long	B	B	18.5"
17.5"	20"	4"	Short	C	A	15.75"
18"	20"	5"	Short	C	A	15.75"
18.5"	22"	5"	Long	B	A	16.75"
19"	20"	7"	Short	C	B	17"
19.5"	22"	7"	Long	B	B	18"
20"	22"	7"	Long	C	B	18"

Standard Features:

- Seat Width 10"-16"
- Seat Depth 11"-16"
- Seat to Floor Height Front 17.50"-20"
- Seat to Floor Height Rear 15.75"-18.5"
- Rear Wheel Options 12", 20" or 22"
- Handrim Option Stainless Steel
- Caster Options 4", 5", 7" or 8"
- Adjustable Height Full Length T-Style Arm (vinyl pad)
- Super Chair Kids Cushion
- Super Chair Kids Back
- Angle Adjustable 10° Bend
- Cane w/ Push bar - High
- Wheel Locks (Push or Pull to Lock)
- Heel Loops
- Aluminum Flip-Up Footplate
- Side Storage Pocket
- Padded Headrest (Kids or Small)
- Choice of Seven Colours (Black, Blue, Burgundy, Orange, Yellow, Red or Silver)
- Round Heavy duty Cross brace Tubing
- Rigid Tilting Frame w/Swingaway Front End

6 REAR WHEEL RIM

Mag Wheels 104RW09 No Charge
 Spoke Wheels 104RW10 \$75

7 REAR WHEEL TIRE

Airless (Solid Rubber) 104RT1 Standard
 Pneumatic 104RT2 \$60

8 HANDRIM OPTIONS

Chrome 104HR1 Standard
 Plastic Coated 104HR2 \$160
 Omit Handrims 105HR3 Less Net \$15

9 ARM TYPE


T-Style Adj. Height Pin Release - Youth 104AR05 Standard
 T-Style Adj. Height Pin Release - Adult 104AR06 \$300
 Cantilever Arm Adj. Angle 104AR11 \$225
 Cantilever Arm Adj. Angle with Guard 104AR12 \$325

Armrest Matrix for -T Style Pin Release Arm

Armrest Style	Transfer Tube Length D1	Arm Guard Length D2	H Range
Adult	6.75	7.50	9.50" - 14.50"
Youth	4.50	5.0	8.75" - 13.75"

Arm Pad Length L

Desk	Full	Curved
10.25"	14.0"	17.0"



Superchair Front STF Height Availability

Rear Wheel	4" Caster					
20"	17.5"	18"				
Rear Wheel	5" Caster					
12"		18.5"				
22"		18.5"				
Rear Wheel	7" Caster					
12"			19"	19.5"		
22"			19"	19.5"		
Rear Wheel	8" Caster					
12"					20"	
22"					19.75"	

10 ARM PAD

<input type="checkbox"/> Desk length Vinyl 10" (pair)	104AP1	No Charge
<input type="checkbox"/> Full length Vinyl 14" (pair)	104AP2	No Charge
<input type="checkbox"/> Contour Gel 17" (pair)	104AP3	Standard

11 BACK OPTIONS

<input type="checkbox"/> Super Chair Kids Back Stock	104BK75	Standard
<input type="checkbox"/> Super Chair Kids Back Non Stock	104BK76	No Charge
<input type="checkbox"/> Super Chair Kids Back Custom	104BK77	\$825
<input type="checkbox"/> Slip on Padded Back Upholstery	104BK48	No Charge
<input type="checkbox"/> Omit Back	104BK49	No Charge

W () x H ()

For other Backs, please fill in the required Backs Order Form and attach it to this order

12 BACK PINS (PRISM BACKS ONLY)

<input type="checkbox"/> Longer Pins for Back Hardware 1" Growth	104BP1	\$50
<input type="checkbox"/> Longer Pins for Back Hardware 2" Growth	104BP2	\$50

13 SEAT OPTIONS

<input type="checkbox"/> Seat Pan	104SE1	Standard
<input type="checkbox"/> Super Chair Kids Cushion Stock	CSK	Standard
<input type="checkbox"/> Super Chair Kids Cushion Non Stock	CNK	No Charge
<input type="checkbox"/> Super Chair Kids Cushion Custom	CCMK	\$685
<input type="checkbox"/> Prism Ideal Cushion Stock	104SC01	No Charge
<input type="checkbox"/> Prism Ideal Cushion Non Stock	104SC02	No Charge
<input type="checkbox"/> Prism Ideal Cushion Custom	104SC03	No Charge
<input type="checkbox"/> Omit Back	104SC20	No Charge

W () x D ()

For other Backs, please fill in the required Backs Order Form and attach it to this order

14 FRONT RIGGING

^A <input type="checkbox"/> 70° Swingaway Upgraded	104FR04	No Charge
^B <input type="checkbox"/> 80° Swingaway Upgraded	104FR19	\$215
^C <input type="checkbox"/> 90° Swingaway Upgraded	104FR05	\$230
<input type="checkbox"/> Elevating Legrest with Swingaway Upgraded	104FR10	\$550
<input type="checkbox"/> Omit Front Rigging	104FR11	Less Net \$20

^A Composite Flip-Up Footplates available for chairs 14" to 16"^B Aluminum Flip-Up Footplates available for chairs 12" to 16"^C Aluminum Angle Adj. Flip-Up Footplates available for chairs 15" to 16"**15 FOOTPLATE TYPE**

<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up <small>Small</small>	104FP13	Standard
<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up <small>Medium</small>	105FP14	No Charge
<input type="checkbox"/> Super Chair Kids Aluminum Flip-Up <small>Large</small>	104FP15	No Charge
<input type="checkbox"/> Composite Flip-Up <small>Small</small>	104FP01	No Charge
<input type="checkbox"/> Composite Flip-Up <small>Medium</small>	104FP06	No Charge
<input type="checkbox"/> Aluminum Adj. Angle Flip-Up <small>Small</small>	104FP10	\$120

16 FOOTREST/LEGREST ACCESSORIES

<input type="checkbox"/> Padded Calf Strap	104FA03	\$70
<input type="checkbox"/> Foam Pad on Footplate (pair)	104FA07	\$70
<input type="checkbox"/> Heel Loops	104FA01	Standard
<input type="checkbox"/> Omit Heel Loops	104FA02	Less Net \$7

17 WHEEL LOCKS

<input type="checkbox"/> Push-to-Lock	104WL1	Standard
<input type="checkbox"/> Pull-to-Lock	104WL2	No Charge
<input type="checkbox"/> Wheel Lock Extension	104WL3	No Charge

18 HEADREST HARDWARE OPTIONS

<input type="checkbox"/> Fixed Hardware	104HE01	Standard
<input type="checkbox"/> Fixed Hardware (+ 4")	104HE02	\$20
<input type="checkbox"/> Single Multi Axial	104HE03	\$158
<input type="checkbox"/> Single Multi Axial (+ 4")	104HE04	\$183
^{D,C} <input type="checkbox"/> Dual Multi Axial	104HE05	\$195
<input type="checkbox"/> Dual Multi Axial (+ 4")	104HE06	\$220
<input type="checkbox"/> Links Multi Adjustable (Comes with 2 links)	104HE10	\$200
<input type="checkbox"/> Links Multi Adjustable Curved (Comes with 2 links)	104HE11	\$240
<input type="checkbox"/> Additional Link	104HE12	\$50
<input type="checkbox"/> 2 Additional Links	104HE13	\$100
<input type="checkbox"/> 3 Additional Links	104HE14	\$150

^D Kids Headrest Foam Pad available for chairs 10"W to 14"W^C Small Headrest Foam Pad available for chairs 15"W & 16"W**19 HEADREST ADAPTOR BRACKETS**

<input type="checkbox"/> Universal Bracket with Clamp (UHRBA)	104HEBA	\$86
---	---------	------

20 HEADREST FOAM PAD OPTIONS

<input type="checkbox"/> Kids Headrest Foam Pad	104HF5	Standard
<input type="checkbox"/> Small Headrest Foam Pad	104HF1	No Charge
<input type="checkbox"/> Small Headrest with Wings Foam Pad	104HF2	\$90
<input type="checkbox"/> Medium Headrest Foam Pad	104HF7	No Charge
<input type="checkbox"/> Medium Headrest with Wings Foam Pad	104HF8	\$90
<input type="checkbox"/> Large Headrest Foam Pad	104HF3	No Charge
<input type="checkbox"/> Large Headrest with Wings Foam Pad	104HF4	\$95
<input type="checkbox"/> Omit Headrest	104HE09	Less Net \$10

21 ANTI TIP

<input type="checkbox"/> Fixed Anti-Tippers	104MA22	No Charge
---	---------	-----------

22 MISCELLANEOUS ACCESSORIES

<input type="checkbox"/> Auto Seat Positioning Strap	104MA03	\$52
<input type="checkbox"/> Velcro Release Positioning Strap	104MA05	\$48
<input type="checkbox"/> 4 Point Harness	104MA56	\$95
<input type="checkbox"/> Quick Release Axle	104AL1	\$110
<input type="checkbox"/> Front and Rear Transit Hook Assembly	104MA30	Standard
<input type="checkbox"/> Padded Lap Tray	104MA09	\$285
<input type="checkbox"/> Half Padded Lap Tray, left side	104MA14	\$350
<input type="checkbox"/> Half Padded Lap Tray, right side	104MA15	\$350

NOTES:**WARRANTY:**

Future Mobility Healthcare warrants the Superchair for kids wheelchair to be free from defects in materials and workmanship for a period of one (1) year from the date of purchase. The frame is warranted for the lifetime of the original purchaser and/or user.