



FUTURE MOBILITY

SEATING & POSITIONING SOLUTION

Manufacturer: Future Mobility

Future Mobility Overbed Table

Future Mobility overbed table are practical, easy-to-use units suitable for hospital and home use. Perhaps you prefer to work with your laptop or simply enjoy your breakfast while in bed. If so, the Medical Overbed Table is the perfect solution. The Future Mobility Overbed Table rises and lowers easily, requiring the slightest of upward pressure, making it user friendly for all ages and needs. Once your desired height is reached, the table top locks firmly and safely into place. The swivel casters make it extremely lightweight and comfortable to move from one position or room to the next. The laminated surface is textured, making it difficult for items to slide off your table. You'll fully appreciate the enormous support and utility this table offers you. If you seek durability and reliability, the Overbed Table is the ideal solution you've been searching for.



The Future Of Mobility - Working For You

"Canadian Made, Assembled, and Owned." In 2000 Future Mobility began designing, manufacturing, and then distributing wheelchairs for the Canadian Home Healthcare market. Our Founder's mission is to help the world by designing high quality, durable mobility products for the North American marketplace.

SPECS

BASE	19.5"W x 28"L
COLOUR	Base is white, support column is powder-coated, table top is a rigid thermofoil material
HEIGHT ADJUSTMENT	Approx. 30" - 45"
PRODUCT WEIGHT	30 lbs
PRODUCT WEIGHT CAPACITY	40 lbs.
WARRANTY	1 year limited

FEATURES

- Light upward touch required for lever and activating the spring-loaded mechanism allowing the table to elevate freely
- Unique H-base design permits use with chairs
- Our hospital grade rigid thermofoil top provides durability and the seamless spill guard edge to reduce infection control issues.
- Durable, heavy duty bases, painted welded tubular steel construction
- Supports up to 40 pounds of evenly distributed weight

CANADA

6750 Professional Court,
Mississauga, Ontario, L4V 1X6

1-888-737-4011

www.futuremobility.ca



**FUTURE
MOBILITY**

USA

1 Buffalo River Place,
Buffalo, NY, 14210

716-783-9130

www.futuremobility.com

