


Serene Prime

THE PRIME CUSHION is designed for users that require mild postural support and offers gentle contours, helping to promote correct seating posture. A seat cushion made of medium density, high resilient, polyurethane foam, that is lightweight, comfortable and with good pressure management.



Stock Sizes – Same Day Shipment

ADP CODE: SE0001225

Part Number	Size	Part Number	Size
PR1616	16" x 16"	PR1818	18" x 18"
PR1618	16" x 18"	PR2016	20" x 16"
PR1816	18" x 16"	PR2018	20" x 18"

Non-Stock Sizes – Shipment within 3 business days

ADP CODE: SE0001225

Part Number	Size	Part Number	Size
PR1620	16" x 20"	PR2220	22" x 20"
PR1820	18" x 20"	PR2416	24" x 16"
PR2020	20" x 20"	PR2418	24" x 18"
PR2216	22" x 16"	PR2420	24" x 20"
PR2218	22" x 18"		



Product Highlights

- Standard water proof cover
- Radiused front edge for foot propelling
- Light weight cushion
- Non slip bottom on washable, incontinent cover
- Promotes mild postural support
- 4 way stretch, water resistant inner liner
- Medium density, high resilient, polyurethane foam

Product Weight: 1.5 Lbs / 0.68 Kg (16"W x20"L)

General Use Cushion

Pressure Reading Scale in mmHG




Pressure Map	Serene Prime
Weight	250 Lbs/113 Kg
Duration	1 hr
Average Pressure	36 mmHg



Polyurethane foam, anatomically shaped



Non slip bottom on washable, incontinent cover



FUTURE MOBILITY

CANADA

6750 Professional Court, Mississauga, L4V1X6

1-888-737-4011

www.futuremobility.ca

USA

11126 Greenstone Ave, Santa Fe Springs, CA 90670

716-783-9130

www.futuremobility.com