

WHEELCHAIR ACCESSORIES DESK LENGTH VINYL ARMREST

FEATURES TEXTURED VINYL FOR A PREMIUM LOOK AND FEEL

THICK HIGH-DENSITY FOAM MEASURING JUST OVER 1 5/8" THICK AT THE CENTER FOR MAXIMUM DURABILITY AND COMFORT

QUALITY PLASTIC BASE WITH METAL THREADED SCREW SLEEVES

LATEX FREE

DESIGN FITS MOST WHEELCHAIRS WITH DESK LENGTH ARMREST PADS

HARDWARE INCLUDED

BLACK VINYL ON A BLACK BASE



UNIQUE DESIGN



LIGHTWEIGHT



TWO (2) YEARS WARRANTY

DESK LENGTH VINYL ARMREST

Wheelchair Armrest Pad in Vinyl – Desk Length

The (desk length) Wheelchair Armrest Pad in Vinyl features premium textured vinyl and super thick high-density foam for durability and all-day comfort.

GET IN TOUCH

orderdesk@futuremobility.ca
www.futuremobility.ca

Tel: 1 888 737 4011
Fax: 1 866 670 3377

6750 Professional Ct, Mississauga,
ON L4V1X6, Canada

